

IMPORTANT NOTIFICATION PRODUCT CHANGE

January, 2015

Dear STERIKING Partner,

Wipak will start implementing the suggested **changes also in the product range of Seal-Sealable (SS) pouches** as well.

1) Change of the indicators

- Current solvent based indicators will be replaced by 3P series of indicators which are water based and lead free complying with ISO 11140-1:2005 class 1.
- The new indicators has been internally validated by Wipak for its safety and performance and also tested and approved by final users in all continents.
- **Color change of indicators will be different** but giving improved contrast and readout; see indicative color chart below.



2) Change of the product print design

- Last year, on the top of feedback and ideas from our customers and partners, we invited the designer resources by www.99designs.com from all corners of the world to propose us a new design on our Steriking packs. The winner was chosen out of 17 competitive works by our global partners and sales teams. Our in-house professionals have now finalized their work and the new design is ready to come out.
- The attached product leaflet informs more details about the new modern outlook.


Start planning now...

... and activate your product change procedure if needed. The changes shall apply to the complete product range of SS- products from the **beginning of February 2015** batch by batch. Product names, codes and packing units will remain unchanged.

If you have any questions or comments, please feel free to contact us.

Sincerely,

WIPAK


 Jaana Kilpinen
 Product Manager
 Steriking, Wipak Health