

FOR IMMEDIATE RELEASE:

CONTACTS:
Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

Clamp for FIS-007

Fraser, MI, March 3, 2022 - Healthmark Industries is delighted to introduce the Clamp for FIS-007 to their $ProSys^{TM}$ Optical Inspection product line. It is an ideal solution for healthcare facilities faced with space constraints. The clamp is designed for mounting the FIS-007 control box to either posts that are 19-38mm in diameter, or rails that are $10 \times 25mm$.

Manufactured from die cast aluminum and stainless-steel, the 130.3 L x 20.8 W x 63.5 Dia. mm clamp can support up to 30 lbs. and has a Clear E Coating finish, which can withstand hospital cleaning solutions. Easily mounted with the use of a Phillips Head screwdriver.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries' mission is to continue to innovate, continue to support, and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.