

STERIKING

DOUBLE PACKING POUCHES

REFERENCE GUIDE FOR MATCHING THE CORRECT INNER SIZE TO OUTER SIZE

IMPORTANT – THIS CHART ONLY OFFERS SUGGESTED PAIRINGS. IT IS CRITICAL TO CONSIDER THE SIZE/BULKLY NATURE OF THE ITEM TO BE PACKAGED

Pouch Item #	Size In.	Fits In	Pouch Item #	Size In.
24	(3 X 6)	➔	5	(4 x 8)
2	(3 x 8)	➔	5	(4 x 10.5)
4	(3 x 10.5)	➔	Roll-42 x 425	(4 x 14)
25	(4 x 6)	➔	6	(6 x 10.5)
3	(4 x 8)	➔	20	(6 x 10.5)
5	(4 x 10.5)	➔	13	(6 x 13.75)
8	(4 x 12)	➔	13	(6 x 15.75)
19	(4 x 14)	➔	28	(6 x 20.5)
12	(4 x 22.5)	➔	R 42 Roll 6" x 600'	
27	(5 x 16)	➔	28	(6 x 20.5)
15	(6 x 8)	➔	10	(8 x 10.5)
6	(6 x 10.5)	➔	10	(8 x 14)
28	(6 x 14)	➔	R43 x 520	(8 x 15.75)
28	(6 x 20.5)	➔	R 43 Roll 8" x 570'	
33	(7 x 12.5)	➔	10	(8 x 15.75)
7	(8 x 10.5)	➔	16	(10 x 15)
10	(8 x 16)	➔	14	(12 x 18)
16	(10 x 15)	➔	21	(12 x 18)'



1-800-521-6224

Samples available for double pouching