



**Instructions for Use: Easy Stringer™**

<b>Brand Name of Product</b>	Easy Stringer™
<b>Generic Name of Product</b>	Easy Stringer™
<b>Product Code Number(s)</b>	ES-100, ES-109, ES-112
<b>Intended Use</b>	To hold the stringer upright in the slotted base to easily load instruments onto the stringer.
<b>Range of Applications for Product</b>	To hold stringers.
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• Solid Plexiglas base to provide weight.</li> <li>• Available sizes: 6, 9- or 12-inch widths</li> <li>• Square base with slot</li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Contents</b>	
<b>Shelf Life</b>	

<b>Instructions for Using Product</b>	
<b>Description of Use(s)</b>	To hold the stringer upright to easily load instruments onto the stringer.
<b>Preparation for Use</b>	
<b>Diagrams (drawings, pictures)</b>	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Place the stringer in the base.</li> <li>2. String one instrument upon another while stacking them on the stringer.</li> <li>3. When the stringer is full, lift out of the base and place in instrument tray.</li> </ol>
<b>Interpretation of Test Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	
<b>Disposal</b>	

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	Clean the Easy Stringer™ with a neutral pH detergent
<b>Cleaning – Automated</b>	
<b>Disinfection</b>	Use Isopropyl alcohol wipe in accordance with the labeling instructions for that wipe.
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	
<b>Storage</b>	
<b>Additional Information</b>	
<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026

1-586-774-7600  
[healthmark@hmark.com](mailto:healthmark@hmark.com)  
hmark.com