


Instructions for Use: LTA Disposable External Brush-Rigid Fibers

Brand Name of Product	LTA Disposable External Brush
Generic Name of Product	Disposable Brush
Product Code Number(s)	SU-510, SU-520
Intended Use	Designed for preliminary cleaning of hard surface items prior to further processing.
Range of Applications for Product	Products that require precision cleaning of devices with hard external surfaces.
Key Specifications of Product	<ul style="list-style-type: none"> Single-use brushes are designed with latex-free thermoplastic handles and nylon bristles. Medium fiber style brush is colored yellow with white bristles and is recommended for routine cleaning. Rigid fiber style brush is colored green with white bristles and is recommended for intensive cleaning.

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	
Packaging Contents	
Shelf Life	

Instructions for Using Product	
Description of Use(s)	Designed for preliminary cleaning of hard surface items.
Preparation for Use	Before each use, check the integrity and the condition of the brush. Replace any worn or damaged brush.
Diagrams (drawings, pictures)	
Steps for Use of Product	<ol style="list-style-type: none"> Consult the manufacturer's IFU for the item you are cleaning and follow those instruction closely. Pay close attention to recommendations on the material composition of the brush trim. Visually inspect box lock area and serrations for any debris / residual / bioburden / blood / protein.
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	<ul style="list-style-type: none"> External brushes are intended for use as cleaning tools only. They are not for use in patient care or aseptic techniques. Wear proper PPE in compliance with facility policy.
Disposal	Because these brushes are used to clean patient used items, it is recommended to dispose of them in a biohazard container.

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	
Disassembly Instructions	
Cleaning – Manual	

Cleaning – Automated	
Disinfection	
Drying	
Maintenance, Inspection and Testing	
Reassembly Instructions	
Packaging	
Sterilization	
Storage	
Additional Information	Intended for single-use. A “use” should be defined by facility policy.
Related Healthmark Products	Brush Cleaning Tools
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com