

Leak Test Tester Log Sheet

Date	Time	Tech's Initials	Department	Olympus Tester Unit MU-1 SN# & PSI 4.0 psi or less is a Fail	ALT-Pro (MAJ-2009) Air Tube #1/ Air Tube #2 SN# & PSI Acceptable Range 2.9 psi-4.1 psi	Olympus Coiled Cable Tester M-155 SN# & PSI Acceptable Range 2.9 psi-4.1 psi	Olympus Handheld SN# & PSI 3.6 psi or less is a Fail	Karl Storz Handheld Tester SN# & PSI 3.0 psi or less is a Fail	Pentax Handheld Tester SN# & PSI 5.0 psi or less is a Fail	Comments/ Actions