

Dear Customer;

In response to your request;

Wipak Medical's sterilization medical grade paper bags are manufactured with medical grade paper, in compliance with EN868-3:1999. We use water resistant and non-corrosive colored adhesives for the bottom seal, and longitudinal seam. Some sizes include a sterilization process indicator and heat seal strip. These meet and comply with the requirements of EN: 868-4:1999.

The medical grade paper bag can be sterilized within wrapped trays and rigid sterilization containers in a standard hospital cycle (steam pre-vacuum; temperature range 270-275°F; duration 4 minutes).

The sterilization parameters assume that all instruments or other items placed within the medical grade paper bag have been cleaned prior to packaging and the process was verified. User must consult with sterilizer manufacture and follow instructions for proper use, cycle parameters, appropriate applications, load capacity, configuration and performance monitoring.

Sincerely,

WIPAK OY



Anne Lehtovuori
Business Manager, Health Care Products
Wipak Medical