



Instructions for Use: Pro-Ties™

Brand Name of Product	Pro-Ties™
Generic Name of Product	Bundling tool
Product Code Number(s)	310-054, 310-053, 310-052, 310-051
Purpose of Product	Pro-Ties™ are used to simplify the process of bundling items. They can significantly reduce time spent sorting.
Range of Applications for Product	
Key Specifications of Product	<ul style="list-style-type: none"> • Stretches to three times its original size • They will not stretch out of shape, rust, or breed bacteria • Steam sterilizable • Reusable • 4 colors: <ul style="list-style-type: none"> ○ 310-054 Yellow ○ 310-053 Blue ○ 310-052 Green ○ 310-051 Red

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	
Packaging Contents	6 ties per box
Shelf Life	

Instructions for Using Product	
Description of Use(s)	Pro-Ties™ are used to simplify the process of bundling items.
Preparation for Use	
Diagrams (drawings, pictures)	
Steps for Use of Product	<ul style="list-style-type: none"> • Place items on Pro-Ties™ • Then lock one end under the other to secure bundle
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	
Disposal	

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	
Disassembly Instructions	
Cleaning – Manual	
Cleaning – Automated	
Disinfection	
Drying	
Maintenance, Inspection, and Testing	

Reassembly Instructions	
Packaging	Can be used within instrument containers, bags, wraps and pouches.
Sterilization	Steam sterilization- 132°C for 4 minutes or 135°C for 3 minutes. Maximum extended time 30 minutes 135°C. Compatible with EO Sterilization.
Storage:	
Additional Information	

Related Healthmark Products	Clamps & Spreaders, Adjustable Stringer, Stringers
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2020-08-05 msmith