

4.8.2015

To whom it may concern

Wipak has implemented a new Steriking 3P design in the beginning of 2015. The new design features for example new sealing line in the end of the pouch. The idea is to demonstrate that the pouch shall be sealed on the top or above of the printed line. This is to avoid that the pouch will be closed from the very end, so that the thumb cut is under the seal. With this line we want to highlight that you should leave some room between the seal and the end of the pouch.

Wipak will make a new 3P leaflet by the end of 2015 to highlight in detail all new features of the Steriking 3P design.

WIPAK OY

Nastola 4.8.2015



Mira Santala
Product Manager/Wipak