



Instructions for Use: Silicone Mesh Covers

Brand Name of Product	Silicone Mesh Covers
Generic Name of Product	Instrument Tray Cover
Product Code Number(s)	SMC-01, SMC-02, SMC-03, SMC-04
Intended Use	To protect delicate instruments by covering instrument trays and restrict movement of items.
Range of Applications for Product	N/A
Key Specifications of Product	<ul style="list-style-type: none"> • Heat resistant up to 135 °C (275 °F). • Flexible frame. • Latex Free. • Silicone cover over stainless steel.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Conditions	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	N/A
Preparation for Use	N/A
Diagrams (drawings, pictures)	N/A
Steps for Use of Product	<ol style="list-style-type: none"> 1. Place instruments in tray. 2. Place Silicone Mesh Cover over the instrument tray. 3. Make sure Silicone Mesh Cover is affixed to restrict instrument movement.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	Clean in accordance with the instrument Mfr.’s IFU. (Note: Silicone Mesh Cover is heat resistant up to 135 °C [275 °F].)
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instruction	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A

Related Healthmark Products	N/A
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026

1-586-774-7600
healthmark@hmark.com
hmark.com