


Instructions for Use: Stainless Steel Tray Tags

Brand Name of Product	Stainless Steel Tray Tags
Generic Name of Product	Tray Tags
Product Code Number(s)	01444, 01452, 01454
Purpose of Product	Alternative for identifying trays, containers, and other key assets.
Range of Applications for Product	N/A
Key Specifications of Product	<ul style="list-style-type: none"> Stainless steel Three styles to choose from

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Con	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	Attach tray tags to trays for identification.
Preparation for Use	N/A
Diagrams (drawings, pictures)	 <div style="display: flex; justify-content: space-around; margin-top: 5px;"> Clip-On Style Round Style Tie-In Rectangular Style </div>
Steps for Use of Product	Three styles for use: <ol style="list-style-type: none"> Clip-on will snap on the rim of the tray. Tie-in round and rectangle tags have pre-cut holes to use ties tags for mounting to trays and other equipment.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	May be processed with the basket of instruments in accordance with the cleaning procedure for those instruments.
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	<ul style="list-style-type: none"> Material compatible with steam sterilization at 132°C (270°F) for 4 minutes, 135°C (275°F) for 3 minutes. Compatible with EO and VH2O2 Sterilization.
Storage	N/A
Additional Information	N/A

Related Healthmark Products	N/A
Other Product Support Documents	Labeling Brochure, Labeling Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com