


Instructions for Use: Stainless-Steel Tray Tags

Brand Name of Product	Stainless-Steel Tray Tags
Generic Name of Product	Tray Tags
Product Code Number(s)	01452, TT-001, TT-002
Purpose of Product	Alternative for identifying trays, containers, and other key assets.
Range of Applications for Product	N/A
Key Specifications of Product	<ul style="list-style-type: none"> Stainless steel Choose from three different styles: <ul style="list-style-type: none"> Clip-on Round Tie-in Rectangular.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Con	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	Attach tray tags to trays for identification.
Preparation for Use	N/A
Diagrams (drawings, pictures)	 <div style="display: flex; justify-content: space-around; margin-top: 5px;"> Clip-on Style Round Style Tie-in Rectangular Style </div>
Steps for Use of Product	<p>Three styles for use:</p> <ol style="list-style-type: none"> Clip-on will snap on the rim of the tray. Tie-in round and rectangle tags have pre-cut holes to use ties/tags for mounting to trays and other equipment.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	May be processed with the basket of instruments in accordance with the cleaning procedure for those instruments.
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	<ul style="list-style-type: none"> Material compatible with steam sterilization:

	<ul style="list-style-type: none"> ○ 132 °C (270 °F) for four (4) minutes ○ 135 °C (275 °F) for three (3) minutes. • Compatible with EO and VH2O2 Sterilization.
Storage	N/A
Additional Information	N/A

Related Healthmark Products	N/A
Other Product Support Documents	Labeling Brochure, Labeling Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com