

Steriking® See-Through Packaging

List of suggested sizes for Double Packing

Inside packaging		Outside packaging	
Article No	Size [mm]	Article No	Size [mm]

Flat Heat-Sealable Pouches

S17	50 x 200	S4	75 x 270
S1	50 x 250	S19	100 x 350
S24	75 x 150	S5	100 x 270
S2	75 x 200	S5	100 x 270
S23	75 x 230	S8	100 x 300
S4	75 x 270	S19	100 x 350
S18	75 x 300	S9	100 x 400
S22	75 x 520	R41	100
S25	100 x 150	S6	150 x 270
S3	100 x 200	S6	150 x 270
S5	100 x 270	S26	150 x 350
S8	100 x 300	S13	150 x 400
S19	100 x 350	S28	150 x 520
S9	100 x 400	S28	150 x 520
S12	100 x 570	R42	150
S27	120 x 400	S28	150 x 520
S15	150 x 200	S2030/S7A	205 x 300/205 x 340
S6	150 x 270	S10	205 x 400
S20	150 x 300	S10	205 x 400
S26	150 x 350	S10A	205 x 440
S13	150 x 400	S11	250 x 500
S28	150 x 520	R43	200
S34	160 x 600	R43	200
S33	180 x 320	S10A	205 x 440
S7	205 x 270	S16	250 x 380
S2030	205 x 300	S30	270 x 440
S7A	205 x 340	S30	270 x 440
S10	205 x 400	S21	300 x 500
S10A	205 x 440	S14	300 x 570
S35	205 x 600	R44/S37	250 / 250 x 700
S2070	205 x 700	R44	250
S16	250 x 380	S21	300 x 500
S11	250 x 500	R45	300
S36	250 x 600	R45	300
S37	250 x 700	R45	300
S29	270 x 350	S21	300 x 500
S30	270 x 440	S14	300 x 570
S21	300 x 500	R46	400
S14	300 x 570	R46	400
S38	320 x 500	R47	400
S31	420 x 500	R49 special	500
S32	420 x 600	R49 special	500

Inside packaging		Outside packaging	
Article No	Size [mm]	Article No	Size [mm]

Flat Heat-Sealable Reels

R39	50	R40	75
R40	75	R41	100
R41	100	R125	125
R125	125	R42	150
R42	150	R43	200
R43	200	R44	250
R44	250	R45	300
R45	300	R46	350
R46	350	R47	400
R47	400	R49 special	500

Gusseted Heat-Sealable Reels

RB50	75	RB51	100
RB51	100	RB51A	120
RB51A	120	RB52	150
RB52	150	RB53	200
RB53	200	RB54	250
RB54	250	RB55	300
RB55	300	RB56	350
RB56	350	RB57	400
RB57	400	RB58 special	450

Self-Sealable Pouches

SS1	90 x 200	SS4	130 x 270
SS2	90 x 250	SS4A	130 x 380
SS4	130 x 270	SS5A	200 x 350
SS4A	130 x 380	SS7	300 x 450
SS5	200 x 350	SS6	250 x 400
SS5A	190 x 330	SS6	250 x 400
SS6	250 x 400	SS7	300 x 450
SS7	300 x 450		

LT-Blueline pouches

LTS7520	75 x 200	LTS1530	150 x 300
LTS1025	100 x 250	LTS1644	160 x 440
LTS1530	150 x 300	LTS2038	205 x 390
LTS2038	205 x 390	LTR44/LTS2550	250 / 250 x 500
LTS2538	250 x 390	LTR45	300
LTS2550	250 x 500	LTR45	300

LT-Blueline Self-Sealable Pouches

LTSS1	90 x 200	LTSS4	130 x 270
LTSS2	90 x 250	LTSS4A	130 x 380
LTSS4	130 x 270	LTSS5A	200 x 350
LTSS4A	130 x 380	LTSS7	300 x 450
LTSS5A	190 x 330	LTSS6	250 x 400
LTSS6	250 x 400	LTSS7	300 x 450
LTSS7	300 x 450		