

Brand Name of Product	Sterile Wipes
Generic Name of Product	Prepackaged Sterile Wipes
Product Code Number(s)	SSW-001, XX6700
Intended Use	For wiping items during/after cleaning or during/after HLD in compliance with the item's IFU.
Range of Applications for Product	For wiping items where sterile, non-linting wipes are indicated in the items IFU.
Key Specifications of Product	<ul style="list-style-type: none"> • Sterile • Absorbent • Non-linting • Non-woven polyester cellulose blend

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	For wiping items where sterile, non-linting wipes are indicated in the items IFU.
Preparation for Use	Follow the item's IFU for wipes use during processing procedures. A partial list of those procedures is below.
Diagrams (drawings, pictures)	
Steps for Use of Product	<p>Possible uses include:</p> <ul style="list-style-type: none"> • Drying the external surfaces after high-level disinfection and before storage. • Covering the distal end and the control section of the endoscope during drying. • Prior to sterilization, drying external surfaces of the item and accessories. • Post disinfection of ultrasound probes in automated disinfection units • Wiping the basin and connectors of an AER to dry. <p>Note: Do not wipe in circular motions. For efficiency, wipe in one direction using parallel, overlapping strokes. Dispose of wipes after each individual use (single-use).</p>
	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	Single-use only.

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A

Sterilization	N/A
Storage	N/A
Related Healthmark Products	Non-linting Wipes
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	ANSI/AAMI/ISO 11137:2013
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com