

Instructions for Use: Emergency Medical Boxes

Brand Name of Product	Emergency Medical Boxes
Generic Name of Product	Medical Boxes
Product Code Number(s)	1310-2 BK, 1410-2 GY, 1610-2 GY, 2010-2 GY, 2310-2 GY, 2710-2, 1611-2 GY, 1711-2 GY, 1711-2 OR, 1960-2 GY, 1910-2 OR, 8040-2 ORANGE
Intended Use	To contain medical supplies for emergencies in a lockable case.
Range of Applications for Product	Storing medical supplies for emergencies.
Key Specifications of Product	<ul style="list-style-type: none"> • Lockable with a breakaway tamper event seal • Pop-up trays / lift-out trays • Various sizes

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	To contain medical supplies for emergencies in a lockable case.
Preparation for Use	Gather medical supplies to place in the Emergency Medical Box.
Diagrams (drawings, pictures)	N/A
Steps for Use of Product	<ol style="list-style-type: none"> 1. Place medical supplies in the Emergency Medical Box. 2. Close the lid and lock the tabs. 3. Insert the appropriate Loks into the provided hole(s) on the box and pull tight.
Interpretation of Test Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	Wipe clean with non-linting cloth using a neutral pH cleaning agent.
Cleaning – Automated	N/A
Disinfection	Use an IPA wipe to disinfect Emergency Medical Box.
Drying	Dry with a non-linting cloth.
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A
Related Healthmark Products	Cynch-Loks, Long Cynch-Loks, Pull Tight-Loks
Other Product Support Documents	Security Products Brochure, Security Products Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600

healthmark@hmark.com
hmark.com