



healthmark
INDUSTRIES CO.
health care products

800-521-6224
www.hmark.com

FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

Healthmark Releases Newest Borescope: FIS-007

August 2, 2021, Fraser, MI - Healthmark Industries is delighted to introduce the FIS-007 to its ProSys™ Optical Inspection product line. It is a borescope designed to visually inspect internal channels of potentially soiled or damaged items with enhanced light, vision, and magnification. It also offers users the option to document with photos and videos the lumens and crevices not visible to the unaided eye.

The FIS-007 features a modular design with interchangeable flexible inspection scope attachments available in diameters of 1.06mm and 1.9mm. These scopes have a working length of 110cm and attach to a 13.3W x 9.9L x 4.7H cm control box, allowing for light level adjustments, image capture*, and video recording*. The 1.06mm scope is designed to inspect internal channels of 1.1mm in diameter or larger, and the 1.9mm scope is designed to inspect internal channels of 2.0mm in diameter or larger.

*These features are only available for the USB control box with FIS-007 software, which is included and installs on Windows 10 PC's

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in delivering surgical instruments and other life-saving medical devices to patients. Healthmark Industries' mission is to continue to innovate, continue to support, and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.