

FOR IMMEDIATE RELEASE:

CONTACTS:
Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

Adapters to Test Light Cord Integrity

Fraser, MI, May 5, 2021 - Healthmark Industries is delighted to introduce the Light Cord Adapters to its ProSysTM Instrument Care product line.

Designed to connect light cords to a light source for inspection, the Light Cord Adapters assist healthcare workers with testing the integrity of light cords for intensity of light and fiber optic damage. Simply, choose the appropriate light cord adapter, attach to the light cord, and secure together. Then, attach a light source to the light cord adapter and secure together, activate the light source and inspect. If any defects or damage are identified, remove the light cord from service.

The Light Cord Adapters are offered for the following makes of light cords: Storz, Stryker, ACMI, and Wolf. The Light Cord Adapters are available for individual purchase. Available light sources are sold separately.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries' mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.