healthmark INDUSTRIES CO. health care products 800-521-6224 www.hmark.com

FOR IMMEDIATE RELEASE:

CONTACTS: Matt Smith Healthmark Industries (586) 774-7600 msmith@hmark.com

Healthmark Offers New UVC LED Box

Fraser, MI October 8, 2020 - Healthmark Industries is delighted to introduce the UVC LED Boxes to its Personal Protection Equipment (PPE) product line. Both boxes contain UV-C LED bulbs that generate ultraviolet beams with wavelengths of 260-280nm, which are effective in sanitizing personal items, such as phones, keys, etc. UVC-001SK has measurements of 216 x 104 x 71mm, and UVC-002SK has measurements of 302 x 269 x 154mm.

For a limited time only, with the purchase of any UVC LED Box, ten Ultraviolet Indicators (UVI-001SK) will be included at no additional cost. These indicators are a great way to verify the function of the UVC LED Boxes. The package of ten comes with a color change reference guide to interpret the test results.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries' mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.