



Instructions for Use: Vinyl Instrument Guards

Brand Name of Product	Vinyl Instrument Guards
Generic Name of Product	Tip Protectors
Product Code Number(s)	09-1-011, 09-1-012, 09-1-013, 09-1-014, 09-1-015, 09-1-016, 09-1-017, 09-1-018, 09-1-019, 09-1-020, 09-2-011, 09-2-012, 09-2-013, 09-2-014, 09-2-015, 09-2-016, 09-2-017, 09-2-018, 09-2-019, 09-2-020, 09-3-031, 09-3-032, 09-3-033, 09-3-034, 09-3-035, 09-03-037, 09-3-040, 09-7-011, 09-7-012, 09-7-013, 09-7-014, 09-7-015, 09-4-031, 09-4-035, 09-4-036, 09-4-038, 09-4-037, 09-4-032, 09-4-033, 09-4-034, 09-6-409, 09-6-129, 09-6-410, 09-6-420
Purpose of Product	Vinyl Instrument Guards are designed to protect items.
Range of Applications for Product	N/A
Key Specifications of Product	<ul style="list-style-type: none"> • Compatible with steam sterilization. • Variety of colors (for Vented and Duo-Vented) • Vented, Duo-Vented, and Cylindrical styles

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Conditions	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	To be placed on the tip of the item.
Preparation for Use	Fit each item with the appropriate size tip.
Diagrams (drawings, pictures)	N/A
Steps for Use of Product	<ol style="list-style-type: none"> 1. Place protector tip on the appropriate size items prior to sterilization. 2. Package items in sterile barrier system as detailed in the item's IFU. 3. Process as indicated in the item's IFU.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	<i>Single-use:</i> Dispose of in compliance with facility guidelines.

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	<p>Approved for sterilization via the following methods:</p> <ul style="list-style-type: none"> • Gravity steam autoclave at 121 °C for 30 minutes. • Pre-vacuum steam autoclave for 4 minutes at 132 °C and 3 minutes at 135 °C. • Maximum steam sterilization conditions: 135 °C for 30 minutes. <p>NOTE: These are minimum cycles.</p>
Storage	N/A

Additional Information	N/A
Related Healthmark Products	N/A
Other Product Support Documents	N/A
Reference Documents	ProSys™ Price List, ProSys™ Brochure
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com