



**Instructions for Use: Wire Mesh & Larger Washer Baskets**

<b>Brand Name of Product</b>	Wire Mesh and Larger Washer Baskets
<b>Generic Name of Product</b>	Washer baskets for medium and larger items
<b>Product Code Number(s)</b>	75.0005.3, 75.0004.6, 75.0005.9, 75.0050.8, 75.0004.8, 75.0004.9, 75.0091.2, 75.0091.1, 72.1400.5, 72.1400.6.
<b>Purpose of Product</b>	To hold parts and items when being washed inside of an automated instrument washer, preventing damage and loss.
<b>Range of Applications for Product</b>	Instruments of from small to larger sizes.
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• Manufactured from stainless steel.</li> <li>• Lid for safe containment of small parts. (Available on the Mesh Washer Basket only)</li> <li>• Designed to allow thorough exposure to cleaning and rinsing solution within an automated instrument washer.</li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	N/A
<b>Packaging Contents</b>	N/A
<b>Shelf Life</b>	N/A

<b>Instructions for Using Product</b>	
<b>Description of Use(s)</b>	Place in the washer basket small items and parts that often get lost or damaged inside of an automated instrument washer.
<b>Preparation for Use</b>	N/A
<b>Diagrams (drawings, pictures)</b>	N/A
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Place items in basket</li> <li>2. Close lid of basket (If applicable)</li> <li>3. Place into instrument washer for cleaning</li> </ol>
<b>Interpretation of Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	N/A
<b>Disposal</b>	N/A

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	N/A
<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	Compatible with steam sterilization. Recommended cycle is 132°C (270°F) for 4 minutes or 135°C (275°F) for 3 minutes in a pre-vac sterilizer.
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A

<b>Related Healthmark Products</b>	Nic-Nac Washer Baskets and Balls
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd.

	Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com
--	---