

## FOR IMMEDIATE RELEASE:

CONTACTS:
Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

## Healthmark Offers New X-Ray Vest and Skirt & X-Ray Apron with Pockets for Cooling Packs

Fraser, MI, November 03, 2021 - Healthmark Industries is delighted to introduce the addition of the X-Ray Vest and Skirt & X-Ray Apron products to its Cools Aids personal cooling products line. Manufactured with a coated nylon fabric exterior and equipped with lead-free core material, the X-Ray Vest and Skirt & X-Ray Apron are designed to keep individuals cool, while also protected from radiation exposure.

The X-Ray Vest and Skirt features both Velcro® and buckles for quick release, with weight being evenly distributed between the shoulders and hips for maximum comfort during long procedures, included are pockets in the front and back of the vest for Cool Pax to help manage the core body temperature of the wearer. The vest has core material that is 0.5 mm lead equivalency at the overlap for front protection and 0.25 mm lead equivalency back protection.

The X-Ray Apron is equipped with Velcro® adjustable shoulder straps to help keep the apron secure. The pockets in the front and back are designed to hold Cool Pax that will help manage the core body temperature of the wearer. The six-inch elastic supports sewn in the back panel helps distribute weight from the shoulders and hips for added back support and comfort. The core material is 0.5 mm lead equivalency for frontal protection.

## About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries' mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit <a href="https://www.hmark.com">www.hmark.com</a> for more information.