


Bougie wire-guided dilators are used primarily in the treatment of esophageal strictures. They are available single-use or reusable. Types include Savary-Gilliard® and American Dilation System dilators.


The dilators are polyvinyl chloride and have a central channel that accommodates a guidewire. They come in various diameters, but the guidewire channel remains consistently the same size. This channel requires brushing and flushing during manual cleaning.

Healthmark's CC2-100-200 Channel Brush is available for this purpose. The nylon bristles and sturdy thermoplastic shaft allows for smooth and effective passage through this lumen.

CC2-100-200 Key Features:

- Single-Use
- Bristle Diameter: 2 mm
- Working Length: 108 cm
- Thermoplastic Shaft
- 12 Individually Packed Per Carton
- Can be used as a push through or pull through brush

 hmark.com

 800.521.6224

 healthmark@hmark.com

