



Healthmark Industries Company

# Personal Protection Equipment

*A complete line of trays, carts, attire and accessories  
for the safe handling, retrieval, decontamination and  
sterilization of instruments and equipment*



SST Tray Systems  
Page 2



SST Tray Systems  
Page 3



SST Accessories  
Page 4



SST Accessories  
Page 4



Transport Carts  
Page 4



Gloves, Gowns & other PPE  
Page 5



Compression Socks  
Page 6



Cool Aids & Scrub Hats  
Page 7



Anti-Fatigue Mat & Humipak  
Page 8

Call:800-521-6224

Fax:586-491-2113

[www.HMARK.com](http://www.HMARK.com)



# SST TRAY SYSTEMS



## SST-2136...The Full Size System

Sized for complete surgical sets. Ideal for retrieval of O.R., O.B. instruments and high volume E.R.'s.

Code No.	Description	Length	Width	Depth	Temp
SST-2136	Amber	21.000	13.000	6.000	285°F
SST-2136 RD	Red	21.000	13.000	6.000	250°F
SST-2136 GKT	Gasket Cover	21.000	13.000	6.000	250°F
SST-2136-SS	Stainless	21.000	13.000	6.000	285°F
SST-2136 LTCH	W/Latch	21.000	13.000	6.000	285°F
SST-2136 RD LTCH	Red w/Latch	21.000	13.000	6.000	250°F



## SST-283...The Mid-Size System

This SST is sized for minor procedure sets including for retrieval of small O.R. and O.B. sets. May be the right size for E.R., Nursing Floors and Clinics. Steristrainer is autoclavable to 285°F. Latches are optional.

Code No.	Description	Length	Width	Depth	Temp
SST-283 RD	Red	18.000	12.000	3.000	275°F
SST-283 RD LTCH	Red w/Latch	18.000	12.000	3.000	275°F



## SST-105...The Small Size System

Our small unit is available with solid or slotted cover, and autoclavable to 275°F. It's ideal for use in Clinics, E.R. and Floors for a small number of instruments.

Code No.	Description	Length	Width	Depth
SST-105 BL	Solid Top-Blue w/Latch	10.000	7.000	6.000
SST-105 RL	Solid Top-RD w/latch	10.000	7.000	6.000
SST-105 RL LTCH	Solid Top-RD w/latch	10.000	7.000	6.000
SST-105 BL LTCH	Solid Top-Blue w/Latch	10.000	7.000	6.000



## Long Saf-T-Soak Systems

These trays are the solution for soaking and cleaning scopes, cameras and other long instruments. Made of chemically resistant plastics, each system includes a base, a cover to limit fumes and a Steri-Strainer to safely remove contents without contacting the solution. Autoclavable to 270°F

Code No.	Description	Length	Width	Depth
SST-2006	Amber Base & Cover	20.875	6.375	8.000
SST-485	Blue Base & Cover	18.000	6.750	6.000
SST-835	Blue Base & Cover	35.500	8.250	6.250

Call:800-521-6224

Fax:586-491-2113

www.HMARK.com

# SST SYSTEMS



## SST-2415 RB

Utilize your own wire sterilization basket as the "steri-strainer" inside this large container. Cover snaps firmly in place. Tray is molded of tough blue polyethylene. Max temp 185°F.

Code No.	Length	Width	Depth
SST-2415 RB	24.000	15.000	5.000



## Disposable SST

Designed for post-procedure transportation of contaminated instruments from procedure areas to the decontamination area, the Disposable SST Tray System is clear and single-use. Manufactured from 100% recyclable plastic (RPET).

Code No.	Length	Width	Depth
DSST-001-10	11.800	12.900	4.000



## Flexible Endoscope Tray

Designed by endoscopy nurses, in this 20" diameter tray flexible scopes can be safely coiled for safe transportation of flexible endoscopes in compliance with industry guidelines. Max temp is 185°F.

Code No.	Length	Width	Depth
2220	20.000		6.000
5-2220	31.150	24.350	
2220 LOK			



## SST-2315

The base tray of this system is large enough to hold a complete instrument tray. The snap down lid and molded-in gasket provide a near watertight seal. Max temp 200°F.

Code No.	Length	Width	Depth
SST-2315	25.000	14.300	6.000



## SST-100

This smallest SST combines convenience with functionality. It has a small footprint but holds a several instruments. The hinged cover snaps shut to contain contents autoclavable to 275 degrees and limit harmful fumes.

Code No.	Length	Width	Depth
SST-100	9.000	4.000	3.000

Call:800-521-6224

Fax:586-491-2113

www.HMARK.com



# SST ACCESSORIES



## SST Tray Stands

With these tray stands, you can replace the outdated basin and ring stand with a system large enough to contain modern instruments - the SST-2136. Choose from a single tray or a double tray style. Manufactured of 100% stainless. Tray systems are sold separately.

Code No.	Length	Width	Height
2136-STND	15.750	28.500	33.000
2136-STND-2	22.000	34.000	36.000

## Vertical Soaker



Ideal to soak M.I.S. type instruments as part of the decontamination process, the Vertical Soaker stands 20" high and 3.5" in diameter. It is graduated to 2000 ml and made of polypropylene.

Code No.	Length	Width	Depth
VS-3520	3.500	0.000	20.000

## 538DP-Applied



Drain contaminated fluids quickly, easily and safely with the optional Drain Plug. Can be installed in any SST Tray. Is autoclavable to 285°F.

Code No.	Description
528DP-AF-AM	Drain Plug - Installed - 1 each, Male & Female Plug System with 5 Foot Hose
538DP	Drain Plug - Installed - 1 each

## Sealed Transport Carts

Covered SST containers in closed carts afford double protection during transporting. Carts are divided into two compartments (except TTCC-20 which is one compartment). Carts allow for delivery of "clean" and pick-up of "dirty" containers at the same time. Constructed of one-piece, double-walled plastic with casters of stainless steel sealed bearings. Completely cart washable. Tray systems are



Code No.	Number & Type of Trays (not included)	Depth	Width	Height
113-20	Small SST Carrier - 4 ea. SST-2136	20.500	27.125	45.000
113-30	Mid Size Cart - 6 ea. SST 2136	20.500	27.125	54.000
113-40	SST Cart - Large - 8 ea. SST-2136	28.500	32.250	45.500
113-50	SST Cart - Extra Large - 12 ea. SST 2136	28.125	32.250	53.250
TTCC-20	SST Cart - Small - 10 ea. SST-283	21.500	32.000	37.500
TTCC-40	SST Cart - Large - 20 ea. SST-283	21.500	33.250	63.750
2220 TRLLY	Flexible Endoscope Transport Trolley	32.150	24.350	50.450



Call:800-521-6224

Fax:586-491-2113

www.HMARK.com





# PPE & Device Covers

## PPE Products

### Disposable Jump Suit

• Non-skid disposable surgical jump suit keeps wearer dry and comfortable. One size fits all and made of impermeable polyethylene for use with sterile gown.

### Disposable Boot Legs

• Made from 2ml Polyethylene that eliminates leakage. An elastic hem holds the bootleg comfortably in place.



Code No.	Description
SHL-001	Face Shield without Drape
SHL-002	Face Shield with Drape
KNEE-HI	Disposable Knee High Leggings
STA-DRI	STA-DRI Disposable Surgical Jump Suit
30125	Jump Suit-White, Elastic Cuff and Wrist
30120	Jump Suit-White Tapered Cuff and Wrist
BOOTLEG	Disposable Bootlegs
42500	Gown with Thumb Loop
VS700CA-S	X-Ray Vest and Skirt -Small
VS700CA-M	X-Ray Vest and Skirt-Med
VS700CA-L	X-Ray Vest and Skirt-Med
VS700CA-XL	X-Ray Vest and Skirt-Med
VS700CA-XXL	X-Ray Vest and Skirt-Med
MN110 SMMD	Decontam Gown-SM-Med
MN110 LGXL	Decontam Gown-Large-XL
MN110 2X3X	Decontam Gown-XXL-XXXL
MN110 4X	Decontam Gown-XXXXL



### Face Shield

• Provides the best splash protection against exposure to bloodborne pathogens.

### Decontamination Gown

• High level strikethrough resistance in the most vulnerable zone, protecting healthcare workers from soil and bioburden.

### X-Ray Vest and Skirt

• The lightweight X-Ray Vest and Skirt & X-Ray Apron comes with a thyroid collar that provides 0.5 mm lead equivalency protection and has an adjustable neck closure with Velcro®.

## Personal Protection Devices

Keep your devices safe and germ-free with our protective covers! We offer disposable, zip lock plastic sleeves designed to keep your phones and tablets clean. We also offer hygienic keyboard covers that protect against pollution or contamination from bacteria and germs. Easy to fit and use, these protective covers may be changed as often as needed to ensure clean use anywhere you go! The protective covers come in two sizes:

(6.25" x 3.38") or (7.88" x 10.25") and keyboard covers fit all standard keyboard sizes.



Code No.	Description	QTY
DPPC-001	Disposable protective phone covers	30
SES-002	Large Phone Disposable Cover	2
DPTC-001	Disposable protective tablet covers	30
2005758	Keyboard Cover Uni Flex Nano	10

Call:800-521-6224

Fax:586-491-2113

www.HMARK.com



# Gloves & Compression Socks



Ultra Long Sleeve Glove



**The Decontam Glove (DSK)** is certified for use in chemotherapy. These gloves are thick enough and tough enough for use in Decontam. This thicker, 15 mil, 12" long textured, powder-free latex glove is the decontam glove of choice by professionals who handle potentially hazardous items.

**The Ultra Long Decontam Glove (ULN16 & ULNCS)** are 16" long powder-free nitrile, latex free available in both 7 and 9ml gloves for safe handling of contaminated items in the Decontam.



**Lined Sleeve Gloves** are a unique fusion of an 11 mil (also available in 15 mil) nitrile glove and a 4 mil protective sleeve. It features 12 inches of sleeve length and an overall length of 28 inches, these gloves provide complete durable hand and arm protection. The sleeve gloves are comfortable, flexible, and lightweight, while the elastic cuff on the sleeve gloves inhibits slippage of the sleeve from the upper arm. These latex-free gloves are a one day use and are resistant to a broad range of hazardous chemicals. Sold in Pairs.

Lined Sleeve Glove

Code No.	Description
DSK-1	Chemotherapy/Decontam, Glove SM - 15 ml
DSK-2	Chemotherapy/Decontam, Glove MED - 15 ml
DSK-3	Chemotherapy/Decontam, LG - 15 ml
DSK-4	Chemotherapy/Decontam, Extra LG - 15 ml
41280	XS Lined Sleeve 15 mil glove - 4 mil Sleeve
41380	Small Lined Sleeve 15 mil glove - 4 mil Sleeve
41480	Med Lined Sleeve 15 mil glove - 4 mil Sleeve
41580	Large Lined Sleeve 15 mil glove - 4 mil Sleeve
41680	X-LG Lined Sleeve 15 mil glove - 4 mil Sleeve

Code No.	Description
ULNCS-1	16 In. 9 Mil Ultra Long Decontam-Sm
ULNCS-2	16 In. 9 Mil Ultra Long Decontam-Med
ULNCS-3	16 In. 9 Mil Ultra Long Decontam-L
ULNCS-4	16 In. 9 Mil Ultra Long Decontam-XL
ULNCS-5	16 In. 9 Mil Ultra Long Decontam-XXL
41270	XS Long Sleeve 11ml latex-free glove/4ml sleeve
41370	Small Long Sleeve 11ml latex-free glove/4ml sleeve
41470	Med Long Sleeve 11ml latex-free glove/4ml sleeve
41570	Large Long Sleeve 11ml latex-free glove/4ml sleeve
41670	X-LG Long Sleeve 11ml latex-free glove/4ml sleeve

## Compression Socks



These comfortable, stylish and customizable unisex Compression Socks are manufactured with moisture wicking, and anti-odor micro nylon fabric.

Designed to help battle muscle fatigue and increase circulation in the lower legs and feet, the Compression Socks reduce swelling and recovery time for the active individual. (Availability of items and sizes subject to change)

Code No.	Description
ZC300ZZ	Below Ankle Mint - S
ZC020ZZ	Below Ankle White - S, M, L
ZC010ZZ	Below Ankle Black - S, M, L
ZN320ZZ	Below Ankle Aqua - S, M
ZC330ZZ	Below Ankle Dark Blue - S
94737-M	Swanky Super Nurse Sock - Med
94699-M	Swanky Healthcare Nurse Sock - Med

## Heart of the Hospital Socks

Manufactured from 80% Acrylic, 20% Nylon and 10% elastic, the Heart of the Hospital Socks are designed to honor CS professionals during Central Service Week. Without the tireless efforts of these individuals, other departments in the hospital wouldn't be able to carry out their necessary lifesaving roles. Made to provide mesh support for the arch and bottom of the foot, the crew length Heart of the Hospital Socks are a great addition to any medical professional's wardrobe.





# Cool Aids, Scrub Hats & Arm Sleeves

## Cool Aids

See more at [StyleYourPPE](#).



Cool Aids are designed to keep your body at a comfortable temperature while you are on your feet for an extended period of time.

You could be in the OR, Endoscopy or the Decontam room for a long shift. When you are wearing Cool Aids under your PPE you will feel more comfortable because the light-weight and the durable fabric is designed for the active individual. These quality items are easy to put on and activate for use. Vests are available with and without cool packs. Neck bands are now available in white with the option to custom print!



Code No.	Description
ECN-001	Cool Aids-Blue Evaporative Cooling Neck Band
ECN-002	Cool Aids-USA Evaporating Cooling Neck Band
ECN-003	Cool Aids-Pink Evaporating Cooling Neck Band
ECN-004	Cool Aids-White Evaporating Cooling Neck Band
ECND-001	Cool Aids-Blue Evaporative Cooling Neck Band-Wide
ECB-001	Cool Aids-Blue Evaporative Cooling Beanie
ECB-002	Cool Aids-USA Evaporative Cooling Beanie
ECB-003	Cool Aids-Pink Evaporative Cooling Beanie
ECS-001	Cool Aids - Blue Evaporative Cooling Skull Cap
ECS-002	Cool Aids - USA Evaporative Cooling Skull Cap
ECS-003	Cool Aids - Pink Evaporative Cooling Skull Cap
VEST-S	Cool Aids-Hybrid Cooling Vest-Small
VEST-M	Cool Aids-Hybrid Cooling Vest-Medium
VEST-L	Cool Aids-Hybrid Cooling Vest-Large
VEST-XL	Cool Aids-Hybrid Cooling Vest-XL
VEST-XXL	Cool Aids-Hybrid Cooling Vest-XXL
VEST-XXXL	Cool Aids-Hybrid Cooling Vest-XXXL
SUCV-001	Single-Use Cooling Vest M/L
SUCV-002	Single-Use Cooling Vest Plus Size
CPX-001	Replacement Cool Packs Inserts (Set of 4)

## Custom Printed Headwear

Disposable Bouffant Caps are manufactured from latex-free polypropylene fabric and have an elastic headband for a secure and comfortable fit. Disposable Scrub Caps are manufactured from spun lace nonwoven rayon material and include a tie-back closure. Both styles are one size fits all and can be custom printed with the design of your choice including your college school colors! Reusable also available.



Code No.	Description
PRSC-WT	Custom Printed Reusable Scrub Hat
PSC-WT	Custom Printed Scrub Hats - Disposable
PBC-WT	Custom Printed Bouffant Caps - Disposable

Call: 800-521-6224

Fax: 586-491-2113

[www.HMARK.com](http://www.HMARK.com)





# ANTI-FATIGUE MAT & HUMIPAK

## Anti-Fatigue Mats

The Anti-Fatigue Mat is designed for healthcare workers in reprocessing areas to help reduce fatigue and discomfort from standing for long hours. Eco-friendly and providing ergonomic comfort, the Anti-Fatigue Mat has a stain and puncture-resistant surface to prevent absorption of dirt and water, a beveled edge at 18 degrees to reduce tripping potential, matte textured top surface to provide optimum friction and anti-microbial properties for added protection. The easy-to-clean slip-resistant Anti-Fatigue Mat is equipped with a high-traction bottom surface certified by the National Floor Safety Institute and is offered in the following sizes: 20 x 32 Inches, 24 x 36 Inches, and 20 x 48 Inches. The Anti-Fatigue Mat is available for individual purchase.

AFM-2032



AFM-2048



AFM-2436



## Keep Soiled Instruments Moist with Humipak

### Time Is the Enemy

Instrument manufacturers, AAMI, AORN and others generally recommend that decontamination of instruments begin within 30 minutes of use so that organic soils, particularly blood, do not dry. But often this is not possible. Is there a solution to extended holding time of soiled instruments?

### Humipak is the Solution

The Humipak consists of a layer of highly absorbent material sandwiched between two layers of water proof film.

### Facilitates Cleaning

To use, place individual instruments, or an entire instrument tray inside the Humipak, add the specified amount of water to the absorbent layer, and seal with the peel away adhesive strip. This creates a water tight, moist atmosphere that will prevent organics from drying over an extended period of time.

### Contents are Visible

The transparent film allows observation of the contents and a list of contents and/instructions can be written directly onto the film.



### Demonstrable Results

Testing has demonstrated that soils will remain moist for up to 3 days, depending upon the type of instruments, and the exact composition of the residue.

### Give Humipak a Try Today

For more information, or to request a free sample to try, call 800-521-6224 or visit [www.hmark.com](http://www.hmark.com) today.

Keep instruments moist and ready for effective cleaning with Humipak!

Code No.	Description	Length	Width	Qty
HPSS2035-25	Small Size	7.900	13.800	PKG-25
HPSS4243-25	Half Din Size	17.000	16.500	PKG-25
HPSS4267	Full Din Size	26.500	16.500	PKG-25
HPSS6577	Double Din Size	30.300	25.500	PKG-10

Call:800-521-6224

Fax:586-491-2113

[www.HMARK.com](http://www.HMARK.com)