




Healthmark Industries Company

ProTech[®] Instrument Trays

A complete line of instrument trays and accessories designed to protect delicate instruments during sterilization, transportation and storage.



Laparoscopy Trays
Page 2



Customizable Trays
Page 4




Instru-Trays
Page 4



Micro ProTech Trays
Page 5


The ProTech[™] System
Delicate instrument damage means lost revenue and expensive repair. The **ProTech[™]** instrument protection system is designed to hold and protect instruments during sterilization, storage and transportation. Choose from the broadest line of trays and combine with plug-in **Secur-It[™]** instrument holders and silicone cushion mats to design a custom tray for each procedure. Trays are made of lightweight, durable materials that are heat and chemically resistant and compatible with steam, EtO and gas-plasma sterilization.



Scope Trays
Page 7



Scope Trays
Page 7



Specialty Trays
Page 10



Instrument Baskets
Page 11

Call: 800-521-6224

Fax: 586-491-2113

www.HMARK.com



Laparoscopic Procedure Trays



Full Size Double Decker Lap Tray

Combine the Laparoscopic Insert with the graduated Secur-It™ holders (3, 5 & 10 mm slots) to accommodate up to 15 lap instruments. Use an optional silicone cushion mat to protect instruments in the bottom - doubling up on available storage space. Or easily install the single width Secur-Its™ to hold additional cannulated instruments. These trays are made of tough, lightweight and corrosion resistant anodized aluminum.

Code#	Length	Width	Height
DDL-1123	23.000	11.000	8.000

Bariatric Laparoscopic Tray


This tray is large enough to accommodate the very long instruments used in bariatric surgeries (max length 65cm). The system is made of tough, lightweight and corrosion resistant anodized aluminum. The laparoscopic Insert with the graduated Secur-It™ holders (3, 5 & 10 mm slots) is designed to accommodate up to 15 lap instruments. The double-decker design holds twice the number of instruments in the same amount of space.

Code#	Length	Width	Height
51527	27.000	11.000	8.000



Cysto Double-Decker Tray

Secur-Its™ can be placed virtually anywhere in tray for maximum efficiency. The double-decker design holds twice the number of instruments in the same amount of space. The graduated Secur-Its™ are designed to hold instruments which are from 3 to 10mm in diameter. Complete tray kit includes 18 single width Secur-Its™. Alternatively, a silicone cushion mat can be substituted for the Secur-Its™ in either the top or bottom tray to protect delicate instruments.



Code#	Length	Width	Height
CYSTO-2214B	20.000	10.000	3.500

Laparoscopic Procedure Trays




Lightweight Poly Double Decker Tray

The two level design of this tray holds double the number of instruments in the same amount of space. This tray includes the lap instrument insert that is outfitted with graduated Secur-It™ instrument holders to accommodate up to 12 laparoscopic instruments. Base tray includes a silicone cushion mat. Tray and components are made of Radel-R®, a medical grade polymer that is lightweight, incredibly strong and compatible with steam, EtO and low temp sterilization.



Code#	Length	Width	Height
4-1130	20.000	10.000	6.000



Extra-Long Lap Tray

This extra long laparoscopic instrument tray holds up to 16 oversized instruments that can be as long as 60cm. The graduated Secur-It™ instrument holders are designed to accommodate instruments from 5 to 10mm in diameter. The included silicone cushion mats protect ancillary items stored in the utility spaces below the laparoscopic instruments. This tray is made of tough but lightweight Radel-R®.

Code#	Length	Width	Height
SDLP-1130	25.500	12.000	5.500

da Vinci® Sterilization Case

Manufactured of tough but lightweight anodized aluminium, this tray is configured to accommodate and secure robotic instrumentation for safe steam sterilization, transport and storage. The Secur-Its™ safely hold and protect up to 8 instruments. The kit includes: 2 ea 10" 15 position universal brackets. Also available in case for 1 or 2 arms.

Code#	Length	Width	Height
31127 DV	23.000	11.000	5.000
300605-DVA	30.000	6.000	2.000






ProTech[®] Trays

Customizeable Full-Size ProTech[™] Trays

Design a ProTech[™] tray to meet your exact needs. Select from several sizes. Options include insert trays, universal instrument brackets, and silicone finger cushion mats. Trays are made of anodized aluminium with mini perforations to allow for proper penetration of sterilant. Universal instrument brackets (refer to price list for details) can be placed virtually anywhere.




Code#	Description	Length	Width	Height
30320-H	2 Level Tray System	20.000	10.000	3.500
31421-H	Universal Tray - Single Layer	21.000	10.500	2.500
31123-T	Universal Tray - Double Layer	22.000	11.000	3.500
31123-T1	Universal Tray - Double Layer	22.000	11.000	5.000
31127-H	Universal Tray - Single Layer	27.000	11.000	3.500



Silicone Finger Mats

Protect delicate surgical instruments with these robust silicone mats. Silicone fingers gently but firmly cushion instruments from damage and help to hold them in place, preventing damage. Choose from a variety of sizes.

Code#	Length	Width	Code#	Length	Width
2-2025	7.750	7.750	2-1055	17.500	7.500
2-1405	7.800	3.400	2-1215	18.300	8.350
SP2313	9.000	5.000	6-1914	19.000	14.000
SP4622	9.000	8.750	2-4005	19.750	9.400
SP2719	10.875	7.625	2-6005	21.100	9.800
2-3215	13.250	8.250	SR10	39.400	9.000



Instru-Trays

- Color-code with five decorator colors.
- Use as Floor Trays, Mayo Trays, Custom Trays, Procedure Trays, Surgical Trays.
- Available with and without holes.
- **New!** Stainless Steel Style.


Code#	Type	Length	Width	Code#	Type	Length	Width
4116H	Plastic	9.000	6.000	HCO-1309	Stainless	13.688	9.813
4106H	Plastic	10.125	6.625	HCO-1510	Stainless	15.125	10.625
41013H	Plastic	13.563	10.438	HCO-1711	Stainless	17.000	11.625
41116H	Plastic	16.000	11.750	HCO-1912	Stainless	19.125	12.625
41317H	Plastic	17.750	13.813	HCO-2116	Stainless	21.250	16.250
41520H	Plastic	20.000	15.000				

Micro ProTech[®] Trays



Metal Micro ProTech[™] Trays

Looking for metal trays that rival plastic trays in pricing, weight and smooth edges? These new 100% metal trays from Healthmark are the answer. Combine these trays with silicon finger mats, or utilize the insert instrument holders to hold instruments in place. All trays are compatible with steam, EtO and Dry-Heat sterilization.



Code#	Description	Length	Width	Height
184014	Base & Cover	6.000	2.500	1.200
184016	Base & Cover	6.500	4.000	1.200
184019	Base & Cover	7.500	4.000	1.200
184021	Base & Cover	7.500	4.000	2.000
181453	Solid (no holes) Tray w/ Holder	11.000	7.000	1.500
184453	Tray w/ Holder	11.000	7.000	1.500
181463	Base, Cover & Holder	7.000	5.000	1.250
182500	Stainless Steel - 5 Clip	7.000	3.333	1.000
182062	Stainless Steel Cassette	7.000	5.500	1.000
182068	Stainless Steel Cassette	11.000	7.000	1.000



Fingermats
and
instrument
holders sold
separately

Hinged-Top Micro Trays

Designed to organize micro-surgical instruments & hardware. Its living hinged cover and snap latch allows fast and secure operation. These trays are designed for use in all forms of sterilization. Use with Silicon cushion mats or plug-in Secur-Its[™] (sold separately) to protect instruments during transportation, sterilization and storage.


Code#	Color	Length	Width	Height
30732-H	Taupe	7.250	3.250	1.125
30842-H	Taupe	8.750	4.500	1.750
31042-H	Taupe	10.250	4.500	1.750
31055-H	White	10.000	5.000	3.250
31462-H	Taupe	14.000	6.250	2.500
31772-H	Blue	17.500	7.750	2.375
50Z900	White	8.000	1.750	1.750



Clear-Top Micro-ProTech[®] Trays

Clear-Top Micro Trays

Organize micro-surgical instruments and accessories in these customizable trays. Removable, clear, snap-on covers make for easy identification of contents. Combine trays with plug-in Secur-Its[™] or silicon finger mats to protect instrument during transportation, sterilization and storage. All trays are compatible with Steam, EtO, and low temperature sterilization. Secur-Its[™] are all designed to fit into the grid system hole pattern. Select from a complete line of Secur-Its[™] inserts with various slot sizes.



Finger mats and Secur-its are optional accessories

The Diminutive De-Lite


Code#	Description	Length	Width	Height
0278	Base & Cover	2.7	1.5	0.8
02781	Base & Cover	2.7	1.5	1.0
027814	Base & Cover	2.7	1.5	1.3



Finger mats and Secur-Its[™] are optional accessories

3 Amigo Micro-ProTechs[™]


Code#	Description	Length	Width	Height
0625	Base & Cover	6.0	2.5	0.8
062514	Base & Cover	6.0	2.5	1.3
6540	Base & Cover	6.5	4.0	0.8
7540	Base & Cover	7.5	4.0	0.8
754015	Base & Cover	7.5	4.0	1.5



Finger mats and Secur-its[™] are optional accessories

Single & Double Decker Twins

Code#	Description	Length	Width	Height
100675	Single Decker Twin	10.0	6.0	0.8
100615	Double Decker Twin	10.0	6.0	1.5
100615B	Insert Tray for 100615	10.0	6.0	0.8




Finger mats are an optional accessory

Two Tier - 2/3 Din Trays

Base and cover with a full tray insert. Manufactured of nearly unbreakable material, these 100% clear trays are ideal for complete sets of instruments that require more than a Micro, but do not need a full-size tray. Combine trays with silicone finger mats to protect instruments during transportation, sterilization and storage.

Code#	Description	Length	Width	Height
101525	Base & Lid	15.0	10.0	1.5
101525B	Full Insert	15.0	10.0	0.8

ProTech® Scope Trays



ER Hand Tray

Code#	Length	Width	Height
30109HN	10.000	9.000	3.000



Key Elevators Box

Code#	Length	Width	Height
30011KE	10.000	5.000	2.000



Currettes Ortho/Neuro Box

Code#	Length	Width	Height
30011CU	10.000	5.000	2.000



Rhoton Suction Box

Code#	Length	Width	Height
30011RS	10.000	5.000	2.000



Bipolar Set Box

Code#	Length	Width	Height
30011BP	10.000	5.000	2.000



Camera Case

Code#	Length	Width	Height
30109CC	10.000	9.000	3.000



Arthroscope Trays

These trays are a tremendous values. The 31042-22 is made of inexpensive polypropylene, while the 3-4110 is made of virtually unbreakable Radel-R®. Steam, EtO, low temp compatible.

Code#	Description	Length	Width	Height
31042-22	Hinged Top Arthro Tray	10.250	4.500	1.750
3-4110	Clear Top Arthro Tray	9.250	3.000	1.330



Osteotome Tray

This case is designed for curved and straight osteotomes, gouges and chisels. It is made of tough but lightweight anodized aluminum. The included Secur-Its™ hold up to 12 instruments from 1/4" to 13/4" diameter.

Code#	Description	Length	Width	Height
SP-3014	Osteotome Tray	9.000	11.000	3.000

Call: 800-521-6224

Fax: 586-491-2113


www.HMARK.com



ProTech[®] Scope Trays



Multi-Scope



Scope Trays with Cord Channel

These innovative trays include Secur-Its™ to safely hold rigid scopes in place and a molded peripheral channel to conveniently hold cables and cords.

Code#	Description	Length	Width	Height
3-5120	Scope Tray w/ Cord Channel	17.000	6.125	2.000
3-5130	Scope Tray w/ Cord Channel	25.000	6.125	2.000

Bariatric Scope Trays

The Bariatric Scope Trays are designed to hold the very long bariatric scopes. Available in two sizes, the wider tray is capable of storing 1-2 scopes and several instruments.

Code#	Description	Length	Width	Height
30126	Bariatric - Single Scope Case	26.000	3.500	3.500
30630	Bariatric - Multi Scope Case	30.000	6.000	2.500



Economy Scope Tray

Designed to safely secure one or two rigid scopes and cord with permanently installed Secur-Its™ designed for 3, 5 & 10mm instruments. Made of tough but economical polypropylene plastic.

Code#	Description	Length	Width	Height
31772-22	Economy Scope Tray	17.500	7.750	2.375

Clear-Top Scope Trays

Lightweight, tough, durable construction, this clear top tray comes with pop-in Secur-Its™ designed to snugly secure 2 rigid scopes..

Code#	Description	Length	Width	Height
261615	Clear Top Scope Tray w/Secure-its	16.000	2.500	1.500
302215	Scope Tray B/L 1 set Secur-Its	22.000	3.000	1.500


ProTech® Scope Trays



Camera Tray

Designed to safely secure expensive cameras and accessories. Manufactured of tough polypropylene plastic, these trays are compatible with EtO, Steam and VH202 sterilization. Included Pop-in Secur-Its™ instrument holders can be easily installed anywhere in the tray for optimal storage. The hinged cover with positive closure means scopes and accessories are well protected.


Code#	Description	Length	Width	Height
31462-22	Camera Tray	14.000	6.250	2.500



Stainless Scope Trays with Secur-Its™

These tough but lightweight scope trays are an excellent alternative to factory supplied trays. Secur-Its™ mount virtually anywhere needed.


Code#	Description	Length	Width	Height
30111-SSC	Scope Tray	11.000	3.000	2.500
30118-SSC	Universal Scope Case	18.000	3.000	2.500



Extra Safe Scope Trays

Secur-Its™ in the base and Secur-Its™ in the cover double protect your scopes - locking them in place with soft silicone.

Code#	Description	Length	Width	Height
3-3131	Two scopes 1mm-10mm	10.670	2.750	1.375
3-3170	One scope 1 mm - 10mm	15.500	2.750	1.375
3-3180	Two scopes 1mm-10mm	15.500	2.750	1.375
3-3200	One scope 1mm - 10 mm	18.000	2.750	1.375



Super Long Scope Tray

Long enough to accommodate the longest rigid scopes. Made of a high performance plastic that is virtually unbreakable, and includes Secur-Its™ to snugly hold 2 scopes.

Code#	Description	Length	Width	Height
9326	Super Long Scope Tray	27.000	4.000	3.000



ProTech[®] Specialty Trays

da Vinci[®] Instrument Tray

The *da Vinci*[®] Instrument Tray is manufactured of anodized aluminium. Its designed for safe transportation and reprocessing of the *da Vinci*[®] instrument.




Tray with 8.5 mm insert



Tray with 12 mm insert

Code#	Description	Width	Length	Height
300605-DVC	DaVinci Instrument Tray	6.00	30.00	5.00

Long Scope Instrument Tray



This Long Scope Instrument Tray is specially designed for long rigid scopes and other delicate instruments.

Manufactured of lightweight durable polypropylene, this tray is compatible with steam, EO and low temperature sterilization. In addition, the inside of the container also features silicone grips that securely hold the scope and instruments in place.

Code#	Description	Width	Length	Height
SMT5510	Container w/lid 1mm-530mm	4.250	26.00	2.50

Fine Mesh Basket with Snap-Closure

Manufactured from stainless steel, the Fine Mesh Basket with Snap-Closure is 105 x 70 x 25 mm and ideal for safely washing small items that often get lost in instrument washers. The small perforations ensure effective water penetration during reprocessing. Simply place items in basket, close lid of basket and place into instrument washer for cleaning. The Fine Mesh Basket with Snap-Closure is available for individual purchase.



Code#	Description	Width	Length	Height
CIL-004	Fine Mesh	10.50	7.00	2.50
CIL-005	Fine Mesh	3.15	1.576	0.788

ProTech® Instrument Baskets



ProTech® Wire Baskets

Choose from several sizes and clever accessories for stainless steel wire baskets for your surgical instruments. They are designed to stack one atop the other for convenient and safe storage of instruments. Basket sizes are DIN-sized.




Code#	Description	Length	Width	Height
75.0005.4	Wire Basket	9.500	9.875	2.000
75.0005.5	Wire Basket	9.500	9.875	2.750
75.0005.6	Wire Basket	18.875	9.875	2.000
75.0004.8	Wire Basket	18.875	9.875	2.750
75.0004.9	Wire Basket	18.875	9.875	4.000
75.0091.2	Cover	9.500	9.875	
75.0091.1	Cover	18.875	9.875	
72.1400.5	Large Washer Basket	23.667	15.750	4.000
72.1400.6	Large Washer Basket	23.667	15.750	7.875



High-Sided Perforated Trays

These trays are made of durable, non-corrosive surgical stainless steel. Designed and manufactured to withstand the demanding environment of surgery and surgical instrument reprocessing.

Code#	Length	Width	Height
HDO-0952	8.875	5.000	2.000
HDO-1062	10.000	6.500	2.000
HDO-1272	12.250	7.750	2.250
HDO-1210	12.250	10.375	2.500
HDO-1610	16.500	10.000	2.500
HDO-2012	20.750	12.750	2.500



Plastic Sterilizing Baskets

These lightweight baskets are considerably less expensive than traditional metal sterilizing trays. Unlike metal, plastic trays have no sharp edges to tear wrap or packaging. Manufactured of durable polypropylene plastic these trays are ideal for the sterilization of all hospital instruments and supplies. Easy to handle and transport. All models are autoclavable to 285°F.


Code#	Length	Width	Height
PIT2015	8.000	6.000	2.000
DINT2525	10.000	10.000	2.000
4106-H	10.125	6.675	1.000
PIT2715	10.640	5.910	1.180
MIT30251	11.800	9.875	2.000
1511-4	15.125	11.250	4.375
DT2-1218	16.000	10.000	3.000
MIT40301	16.500	12.000	3.000
DT2-2136	17.750	10.750	4.000
1511-7	19.750	19.750	3.750
DINT4825	20.000	10.000	2.000



ProTech[®] Baskets & Mats

Perforated Stainless Steel Trays


This superb range of perforated sterilization baskets are made with a strong stainless steel frame and ensure a construction to last a lifetime. The perforated design ensures maximum steam penetration but limits the chance for delicate instrument points to become damaged. DIN sized trays for compatibility with Closed Containers and Automated Instrument Washers.



Code#	Length	Width	Height
HBO-0520	20.00	5.00	3.50
PSST-009	10.00	6.50	2.50
PSST-007	10.50	10.00	3.50
PSST-005	10.50	15.00	3.50
PSST-002	20.00	5.00	3.50
PSST-001	15.50	5.00	3.50
PSST-012	20.25	13.00	3.50
PSST-003	20.25	13.00	3.50
PSST-011	24.00	6.00	3.50
PSST-004	25.00	15.00	3.50
PSST-013	25.25	8.00	3.50
PSST-008	26.00	9.00	3.50
PSST-006	30.00	6.00	3.50

Magnetic Mat

Designed with strong round-raised magnets, the 12" x 16" purple silicone Magnetic Mat helps prevent items from falling, allowing healthcare professionals to grasp for easy removal. Intended to be utilized as a hands-free "neutral-zone" sharps transfer area, the reusable waterproof mat conforms to all surface shapes and can be fastened to drapes with a standard towel clamp for ease of access. The Magnetic Mat is manufactured with cut-lines that can be split in half for two smaller sizes and is heat resistant up to 275°F



Code#	Description	Width	Length
11116	Magnetic Mat	12.00	16.00