



healthmark  
INDUSTRIES CO.  
health care products

800-521-6224  
www.hmark.com

## **FOR IMMEDIATE RELEASE:**

### CONTACTS:

Matt Smith  
Healthmark Industries  
(586) 774-7600  
msmith@hmark.com

### **Healthmark Industries Offers Steriking® Pouches for Long Medical Devices**

Fraser, MI April 27, 2022 -Healthmark Industries is delighted to introduce the Steriking® Pouches for Long Medical Devices to its Steriking® Sterilization product line. Manufactured with heat-sealed upper corners that prevent dust from collecting at the opening seal to reduce the risk of contamination when opened, the Steriking® Pouches for Long Medical Devices are designed as packing material for sterilization of long devices by steam and are offered in the following sizes: 7.875W x 31.5L In. and 10W x 35.5L In.

The single-use Steriking® Pouches for Long Medical Devices are constructed of heavy weight 100 gsm paper that is heat-sealed together with a multi-ply PET/PE-plastic laminate. Equipped with strong and durable multiline parallel seams that promote a clean peel, the see-through heat-sealable pouches allow double pouching and have thumb notches at each end, which make the filling and opening of a pack easier.

### ***About Healthmark Industries Co., Inc.***

*Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other lifesaving medical devices to patients. Healthmark Industries' mission is to continue to innovate, continue to support, and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit [www.hmark.com](http://www.hmark.com) for more information.*