

Tracheostomy Brushes



The tracheal channel brushes are constructed of low absorbent thermoplastic fibers encapsulated within a twisted curved stainless steel frame and attached to a super flexible thermoplastic handle/wand.

The curvature of the brushes allows for confident and easy lumen penetration, while the flexible handle offers maximum maneuverability during the cleaning process.

Key Features:

- Single-Use
- Curved Bristle
- Flexible Handle/Wand
- Overall Length: 15 cm
- 4 Different Diameters Available
- 10 Per Package



Item Number	Brush Diameter	Channel Diameter	Product Image
SUCT01-10	8 mm	6 mm - 7 mm	
SUCT02-10	10 mm	8 mm - 9 mm	
SUCT03-10	11 mm	10 mm	
SUCT04-10	13 mm	11 mm - 12 mm	
SUCTPANA-40	Includes 10 of each of the above items, for a total of 40 Tracheostomy Brushes		