




Brand Name of Product	Hanging Tag
Generic Name of Product	Labels
Product Code Number(s)	04172
Purpose of Product	For identifying trays, containers, and other key assets.
Range of Applications for Product	The surface of the tag is ideal for barcode labels.
Key Specifications of Product	<ul style="list-style-type: none"> • Composition: Paper, 224 g/m² • Color: White • Size: 3.5 x 1.9 x 0.78 in • Non-adhesive and heat resistant up to 135 °C (275 °F)

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	Case of one thousand (1,000)
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	For identifying trays, containers, and other key assets.
Preparation for Use	N/A
Diagrams (drawings, pictures)	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Tag items being identified in the tray or container 2. Secure tag by tying the precut slots to either the tray, container, or other equipment. (Note: For security, use a lock when using closed containers.)
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	Single use
Disposal	N/A

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A

Sterilization	<ul style="list-style-type: none"> • Steam sterilization, up to 135 °C (275 °F) for 3 minutes and 132 °C (270 °F) for 4 minutes, • Materially compatible with EO sterilization.
Storage	N/A
Additional Information	N/A

Related Healthmark Products	N/A
Other Product Support Documents	Label Price List, Label Brochure
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com