

## Instructions for Use: Wall Mount for Interchangeable Magnifying Lamp

<b>Brand Name of Product</b>	Wall Mount for Interchangeable Magnifying Lamp
<b>Generic Name of Product</b>	Wall Mount
<b>Product Code Number(s)</b>	MAG-WM
<b>Intended Use</b>	To hold the Interchangeable Magnifying Lamp on a vertical surface for optical inspection of items.
<b>Range of Applications for Product</b>	Healthmark's Interchangeable Magnifying Lamp.
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• Permanent mounting on vertical surfaces.</li> <li>• B-Wall mount.</li> <li>• Materials:               <ul style="list-style-type: none"> <li>○ Steel with an enamel coating</li> <li>○ Color: White.</li> </ul> </li> </ul>

Shipping & Storage	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	N/A
<b>Packaging Contents</b>	N/A
<b>Shelf Life</b>	N/A

Instructions for Using Product	
<b>Description of Use(s)</b>	Holds an Interchangeable Magnifying Lamp steady for optical inspection of items.
<b>Preparation for Use</b>	Tools needed: <ul style="list-style-type: none"> <li>• Philips screwdriver</li> <li>• Three (3) screws to attach Wall Mount to a vertical surface.</li> </ul>
<b>Diagrams (drawings, pictures)</b>	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Place B-Bracket to a vertical surface. (Note: Consult with your facilities management department to locate the wall studs or installing for fastening brackets.</li> <li>2. Using a Philips screwdriver, place the screws (an appropriate sized screw) in the three (3) holes of the B-Bracket and tighten to secure to the vertical surface.</li> <li>3. Insert the arm of the Interchangeable Magnifying Lamp into the top opening of the B-Bracket.</li> </ol>
<b>Interpretation of Test Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	N/A
<b>Disposal</b>	N/A

Reprocessing Instructions	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	Use with Isopropyl Alcohol wipes per the wipe manufacturer's IFU.

<b>Cleaning – Automated</b>	N/A
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	N/A
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A
<b>Related Healthmark Products</b>	Interchangeable Magnifying Lamp
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>