



<b>Brand Name of Product</b>	Tray Identification Tags
<b>Generic Name of Product</b>	Tray Identification Tags
<b>Product Code Number(s)</b>	RING001, TG001, TG001-E, TG011, TG012, TG013, TG014, TG015, TG016, TG017, TG018, TG019, TG020, TG021, TG022, TG023, TG024, TG025, TG026, TG011-E, TG012-E, TG013-E, TG014-E, TG015-E, TG016-E, TG017-E, TG018-E, TG019-E, TG020-E, TG021-E, TG022-E, TG023-E, TG024-E, TG025-E, TG026-E, TIE003, STT-001,
<b>Intended Use</b>	For identification, tracking, and tracing purposes.
<b>Range of Applications for Product</b>	A range of colored plastic labeling tags for identification/tracking of trays, baskets, and other items.
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• Resin bonded.</li> <li>• Various colors.</li> <li>• X-ray detectable.</li> <li>• Autoclavable.</li> </ul>

Shipping & Storage	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	N/A
<b>Packaging Contents</b>	N/A
<b>Shelf Life</b>	N/A

Instructions for Using Product	
<b>Description of Use(s)</b>	For use with trays, baskets, and other items for identification purposes.
<b>Preparation for Use</b>	N/A
<b>Diagrams (drawings, pictures):</b>	
<b>Steps for Use of Product</b>	<b>HOW TO USE TRAY TAGS</b>

	<p><b>(Note: Care should be taken to ensure the writing (with a suitable marking pen) will not be removed or bleed during processing. We recommend to test this prior to placing tag in use.)</b></p> <p><b>Writing directly onto the Tray Identification Tag:</b></p> <ol style="list-style-type: none"> <li>Information may be written directly on the Tag.</li> <li>Either or both sides of the Tag may be written on. <ol style="list-style-type: none"> <li>Use a pre-printed self-adhesive label.</li> <li>Pre-printed self-adhesive labels may be stuck onto either side of a Tag.</li> </ol> </li> </ol> <p><b>When Attaching the Tray Identification Tag to Tray or Basket or other items:</b></p> <ol style="list-style-type: none"> <li>Use the holes in the Tag to attach it to the tray or basket.</li> <li>Use a Split Ring (RIN001) or a Security Tie (TIE003) for this purpose.</li> <li>If X-ray detectability is required, use a Split Ring (RIN001) to attach the Tag to the tray since the split ring is X-ray detectable.</li> </ol>
<b>Interpretation of Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>Facilities should develop a policy for the placement of the Tray Identification Tag on the tray or basket. This should include a risk assessment to avoid potential interference with effective processing or aseptic presentation.</li> <li> Store away from any sources of ignition.</li> </ul>
<b>Disposal</b>	<ul style="list-style-type: none"> <li>Users should reference facility procedures.</li> <li>Unused Tray Identification Tags may be disposed of in the plastic waste receptacle.</li> <li>Used Tray Identification Tags should be disposed of by considering possible microbial contamination if not decontaminated before disposal.</li> </ul>

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	N/A
<b>Cleaning – Automated</b>	Tray Identification Tags are compatible with washer disinfectors conforming to ISO 15883 series of standards.
<b>Disinfection</b>	N/A
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	<ul style="list-style-type: none"> <li>Materially compatible with Steam and EO.</li> <li>Recommended cycle for steam sterilization: <ul style="list-style-type: none"> <li>132 °C (270 °F) for four (4) minutes</li> <li>135 °C (275 °F) for three (3) minutes.</li> </ul> </li> <li>Extended steam sterilization cycle: 135 °C (275 °F) for 30 minutes.</li> </ul>
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A

<b>Related Healthmark Products</b>	N/A
<b>Other Product Support Documents</b>	Labeling Brochure, Labeling Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com