

Brand Name of Product	UnderGuard™ for VH202
Generic Name of Product	Tray Liners
Product Code Number(s)	MT-2030, MT-1520, MT-1224
Purpose of Product	For use under the tray to protect outer wrap from sharp corners and imperfections of trays.
Range of Applications for Product	N/A
Key Specifications of Product	<ul style="list-style-type: none"> Manufactured from spun polypropylene. Unique purple color helps distinguish protective layer

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	<ul style="list-style-type: none"> Temperature: 15 °C (59 °F) to 30 °C (86 °F) Humidity: 30–70% rh.
Packaging Contents	N/A
Shelf Life	Five (5) years from manufacturer's date.

Instructions for Using Product	
Description of Use(s)	To protect outer wrap from sharp corners and imperfections of trays.
Preparation for Use	<ul style="list-style-type: none"> Select a liner which will best fit the target basket. <i>Wrapped tray</i>- the liner must be inside the wrap and may be used inside or under the basket. <i>Closed container</i>- the liner must be used inside the basket (either inside or below the basket).
Diagrams (drawings, pictures)	
Steps for Use of Product	<ol style="list-style-type: none"> Place a liner either inside or beneath basket. (Note: In the event you are using a container, the liner should be placed in basket first.) Place items in tray. Wrap or close package. Place in VH2O2, VH2O2-O3, or steam sterilizer.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	MT-XXXX are number 5 ♻ for recycling.

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	Material compatible with VH2O2/O3 and steam sterilization up to 137 °C (278.6 °F).
Storage	<ul style="list-style-type: none"> Temperature: 15 °C (59 °F) to 30 °C (86 °F) Humidity: 30–70% rh.
Additional Information	N/A

Related Healthmark Products	N/A
------------------------------------	-----

Other Product Support Documents	Sterilization Brochure, Sterilization Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com