

### Laparoscopic Procedure Trays

Code No	Description	Length	Width	Height	Mat Code No	Length	Width
CYSTO-2214B	Cysto Tray	20.000	10.000	3.500			
51527	Double Decker Bariatric	25.620	10.700	7.660			
300605-DVA	Da Vinci Arm Sterilization Case	30.000	6.000	5.000			

### Customizable Protech® Procedure Trays

Code No	Description	Length	Width	Height	Mat	Length	Width
30320H	2 Level Tray System	20.000	10.000	3.500	2-1215	18.000	8.000
31421-H	Universal Tray - Single Layer	21.000	10.500	2.500	2-4005	19.750	9.400
31123-T	Universal Tray - Double Layer	22.000	11.000	3.500	2-6005	21.100	9.800
31123-T1	Universal Tray - Double Layer	22.000	11.000	5.000	2-6005	21.100	9.800
31127-H	Universal Tray - Single Layer	27.000	11.000	3.500			

### Instrument Brackets for Customizable Protech® Trays

Code No	Description	Length
P-01	Post	
TRI-CUT BRCKT	1" universal bracket 3, 5 & 10 mm	1.000
56-UB	5" 6 position universal bracket	5.000
710-UB	7.5" 10 position universal bracket	7.500
1015-UB	10" 15 position universal bracket	10.000

### Instru-Trays

Code No	Material	Length	Width	Code No. w/o Holes
4116H BL	Plastic	9.000	6.000	4116 BL
4106H	Plastic	10.125	6.625	4106
41013H BL	Plastic	13.563	10.438	41013 BL
41116H BL	Plastic	16.000	11.750	41116 BL
41317H BL	Plastic	17.750	13.813	41317 BL
41520H BL	Plastic	20.000	15.000	41520 BL
HCO-1065	Stainless	10.000	6.500	
HCO-1309	Stainless	13.688	9.813	HCO-1309-NP
HCO-1510	Stainless	15.125	10.625	HCO-1510-NP
HCO-1711	Stainless	17.000	11.625	
HCO-1912	Stainless	19.125	12.625	HCO-1912-NP
HCO-2116	Stainless	21.250	16.250	

### Metal Micro-ProTech® Trays

Code No	Description	Length	Width	Height
184014	Base & Cover	6.000	2.500	1.200
184016	Base & Cover	6.500	4.000	1.200
184019	Base & Cover	8.000	4.000	1.200
184021	Base & Cover	8.000	4.000	2.000
181453	Solid (no holes) Tray w/ Holder	11.000	7.000	1.500
184453	Tray w/ Holder	11.000	7.000	1.500
181463	Base, Cover & Holder	7.000	5.000	1.250
182500	Stainless Steel - 5 Clip	7.000	3.333	1.000
182062	Stainless Steel Cassette	7.000	5.500	1.000
182068	Stainless Steel Cassette	11.000	7.000	1.000

### Hinged-Top Micro-ProTech® Trays

Code No	Description	Length	Width	Height	Mat Code No
30732-H	Taupe	7.250	3.250	1.125	307-M
30842-H	Taupe	8.750	4.500	1.750	308-M
31042-H	Taupe	10.250	4.500	1.750	310-M
31055-H	White	10.000	5.000	3.250	310-M
31462-H	Taupe	14.000	6.250	2.500	314-M
31772-H	Blue	17.500	7.750	2.375	317-M
50Z900	White	8.000	1.750	1.750	

### Specialty Micro-ProTech® Trays

Code No	Description	Length	Width	Height
30109HN	ER Hand Tray	10.000	9.000	3.000
30109CC	Camera Case	10.000	9.000	3.000
30011KE	Key Elevators Box	10.000	5.000	2.000
30011CU	Currettes Box	10.000	5.000	2.000
30011RS	Rhoton Suction Box	10.000	5.000	2.000
30011BP	Bipolar Set Box	10.000	5.000	2.000

### Metal Hinged Top Micro Trays


Code No	Description	Length	Width	Height	Mat	Length	Width
31020-H	Tray & Cover & Insert	20.000	10.000	2.650	2-1215	18.000	8.000

## Clear-Top Micro-ProTech® Trays

Code No	Description	Length	Width	Height	Mat
0278	Base & Cover	2.7	1.5	0.8	0278-M
02781	Base & Cover	2.7	1.5	1.0	0278-M
027814	Base & Cover	2.7	1.5	1.3	0278-M
0625	Base & Cover	6.0	2.5	0.8	0625-M
062514	Base & Cover	6.0	2.5	1.3	0625-M
6540	Base & Cover	6.5	4.0	0.8	6540-M
7540	Base & Cover	7.5	4.0	0.8	7540-M
754015	Base & Cover	7.5	4.0	1.5	7540-M
100675	Single Decker Twin	10.0	6.0	0.8	100675-M
100615	Double Decker Twin	10.0	6.0	1.5	100615-M
100615B	Insert Tray for 100615	10.0	6.0	0.8	100615BM
101525	Base & Lid	15.0	10.0	1.5	101525-M
101525B	Full Insert	15.0	10.0	0.8	101525B-M

## Silicone Plug-In Secur-It® for Micro-ProTech™ Trays

Code No	Description	Height	Length
2-1/4" wide by 5/8" high Pop-In Secur-Its: Fits all sizes of Micro Protech's			
A237591K	One 4mm slot, two plugs	0.625	2.250
A237592K	Two 4mm slots, two plugs	0.625	2.250
A237593K	Three 4mm slots, two plugs	0.625	2.250
A237593	Three 1mm slots, two plugs	0.625	2.250
A237596	Six 1mm slots, two plugs	0.625	2.250
3-3/4" wide by 5/8" high Pop-In Secur-Its: Fits all sizes of Micro Protech's except #30732-H and #0625			
A387596	Six 1mm slots, three plugs	0.625	3.750
A3875912	Twelve 1mm slots, three plugs	0.625	3.750
A387592	Two 4mm slots, three plugs	0.625	3.750
A387593	Three 4mm slots, three plugs	0.625	3.750
A387594	Four 4mm slots, three plugs	0.625	3.750
A387595	Five 4mm slots, three plugs	0.625	3.750
2-1/2" wide by 1-1/8" high Pop-In Secur-Its: Fits all sizes of Micro Protech's except #30732-H, #0625, #7540, #100675, and #100615B			
A2616901	One 12mm, one 8mm, two plugs	1.125	2.500
A2616902	Two 8mm slots, two plugs	1.125	2.500
2-1/2" wide by 1-1/8" high Pop-In Secur-Its: Fits all sizes of Micro Protech's except #30732-H, #0625, #7540, #100675, and #100615B			
A32759P2	One 5mm slot, three plugs	1.125	3.750
A327593	Three 5mm slots, three plugs	1.125	3.750
A327593L	Two 12mm, one 5mm slot, three plugs	1.125	3.750
A32759P1	One 12mm slots, three plugs	1.125	3.750




## ProTech® Scope Trays

Code No	Description	Length	Width	Height
30126	Bariatric - Single Scope Case	26.000	3.500	3.500
30630	Bariatric - Multi Scope Case	30.000	6.000	2.500
30111-SSC	Scope Tray	11.000	3.000	2.500
30118-SSC	Universal Scope Case	18.000	3.000	2.500
31042-22	Hinged Top Arthro Tray	10.250	4.500	1.750
31772-22	Economy Scope Tray	17.500	7.750	2.375
261615	Clear Top Scope Tray	16.000	2.500	1.500
302215	Scope Tray B/L 1 set Secur-Its	22.000	3.000	1.500
31462-22	Camera Tray	14.000	6.250	2.500
300605-DVC	Da Vinci Lens Case & Lid	30.000	6.000	5.000
SMT5510	Container w Lid 1mm - 530mm	26.000	4.250	2.500



## Silicone Finger Mats & Magnetic Mat

Code No	Description	Length	Width
2-2025	Silicone Mat - Bubble Bottom	7.750	7.750
2-1405	Silicone Mat - Bubble Bottom	7.800	3.400
SP2313	Silicone Mat - Perforated Base - Pink	9.000	5.000
SP4622	Silicone Mat - Perforated Base - Pink	9.000	8.750
SP2719	Silicone Mat - Perforated Base - Pink	10.875	7.625
2-3215	Silicone Mat - Bubble Bottom	13.250	8.250
2-1055	Silicone Mat - Bubble Bottom	17.500	7.500
2-1215	Silicone Mat - Bubble Bottom	18.300	8.350
6-1914	Silicone Mat - Bubble Bottom (no holes)	19.000	14.000
2-4005	Silicone Mat - Bubble Bottom	19.750	9.400
2111-M	Silicone Mat - Bubble Bottom	21.000	11.000
2-6005	Silicone Mat - Bubble Bottom (2 piece)	21.100	9.800
SR10	Silicone Mat - Perforated Base - Roll - Pink	39.400	9.000
11116	Magnetic Mat	16.000	12.000



## ProTech® Wire Baskets


Code No	Description	Length	Width	Height
75.0005.4	Wire Basket	9.500	9.875	2.000
75.0005.5	Wire Basket	9.500	9.875	2.750
75.0005.9	Wire Basket	18.875	9.875	2.000
75.0005.6	Wire Basket	18.875	9.875	4.000
75.0091.2	Cover for Wire Basket	9.500	9.875	
75.0194.2	Cover Wire Basket	18.875	9.875	
72.1400.5	Large Washer Basket	23.667	15.750	4.000
72.1400.6	Large Washer Basket	23.667	15.750	7.875
7501542	Secur-Its -12mm Dia.	1.440		3.969
7501543	Secur-Its - 1-10mm Dia.	1.440		3.969
7501544	Secur-Its - 12.5mm, 6.7, 4.3 Dia.	1.440		3.969
7501545	Secur-Its - 5mm, 3mm, 1.8mm Dia.	1.440		3.969
CIL-004	Fine Mesh Basket with Snap-Closure	10.500	7.00	2.500

## High-Sided Perforated Trays

Code No	Description	Length	Width	Height
HDO-0952	Perforated Tray	8.875	5.000	2.000
HDO-1062	Perforated Tray	10.000	6.500	2.000
HDO-1272	Perforated Tray	12.250	7.750	2.250
HDO-1210	Perforated Tray	12.250	10.375	2.500
HDO-1610	Perforated Tray	16.500	10.000	2.500
HDO-2012	Perforated Tray	20.750	12.750	2.500

## Plastic Sterilizing Baskets

Code No	Length	Width	Height	Qty
PIT2015	8.000	6.000	2.000	EACH
DINT2525	10.000	10.000	2.000	EACH
4106H	10.125	6.675	1.000	EACH
PIT2715	10.640	5.910	1.180	EACH
MIT30251	11.800	9.875	2.000	EACH
1511-4	15.125	11.250	4.375	EACH
DT2-1218	16.000	10.000	3.000	EACH
MIT40301	16.500	12.000	3.000	EACH
DT2-2136	17.750	10.750	4.000	EACH
1511-7	19.750	19.750	3.750	CS-6
DINT4825	20.000	10.000	2.000	EACH



## Perforated Stainless Steel Trays

Code No	Length	Width	Height
PSST-009	10.00	6.50	2.50
PSST-007	10.50	10.00	3.50
PSST-005	10.50	15.00	3.50
PSST-010	10.50	20.00	3.50
PSST-001	15.50	5.00	3.50
PSST-002	20.00	5.00	3.50
PSST-012	20.25	13.00	3.50
PSST-003	20.25	13.00	3.50
PSST-011	24.00	6.00	3.50
PSST-004	25.00	15.00	3.50
PSST-013	25.25	8.00	3.50
PSST-008	26.00	9.00	3.50
PSST-006	30.00	6.00	3.50

