



Sterilizer Cleaning Brushes

A curated assortment of brushes designed for cleaning autoclaves and sterilizers

Designed to be used with a suitable cleaning solution, Healthmark's Sterilizer Cleaning Brushes assist with the internal cleaning of autoclaves and sterilizers for which cleaning of surfaces is a recommended step.

hmark.com
800.521.6224

Utility Scrub Brush / USB-01

Ergonomic handle supports a strong yet comfortable grip while brushing

Plastic handle is lightweight & break-resistant

Crimped bristles resist absorbing moisture and wash easily

BPA Free



Double-Ended Acetabular Reamer Brush / 545

Doubles as a drain strainer cleaning brush



Telescoping Handle / TH-01

Handle easily twists into the desired length - Ranges from 54" to 8'

Fiberglass handle includes a keyhole at the top edge for hanging



Deck Brush / DB-01

Telescoping Handle Sold Separately

Foam-injected block provides durability

Side bristles allow for dual cleaning of floor and corners

Crimped bristles resist absorbing moisture and wash easily

Standard 3/4" threaded handle



Tank Brush / TB-01

Telescoping Handle Sold Separately

Made on a chemical-resistant polypropylene block

Bristles flare to give you access to all corners

Crimped bristles resist absorbing moisture and wash easily

Standard 3/4" threaded handle