



Outside-of-Tray Count Sheet Holder Solutions

Healthmark offers a variety of solutions to help your department easily display count sheets on the outside of wrapped trays and closed containers. Options for wrapped trays include the adhesive-backed PBR-46 and PBR-76. Options for closed containers include the Slotted Tray Label, available in white or pink, and the Stainless Steel Count Sheet Holder.

hmark.com

800.521.6224



Slotted Tray Label

White 403225 SLTD / Pink 403225 SLTD PK

The self-looping Slotted Tray Label allows you to attach a rolled up surgical instrument count sheet on the outside of a closed sterilization container.



Stainless-Steel Count Sheet Holder

CSH-001


The Stainless-Steel Count Sheet Holder is designed to secure and hold folded count sheets while clipped on the outside of closed container systems. Simply fold the count sheet, place it securely within the stainless-steel loops, and use the clip portion to attach it to the outside of the container.



PBR-76 Count Sheet Holder

PBR-76

The cellulose PBR-76 Count Sheet Holder is designed to encase the count sheet and be adhered to the outside of a wrapped tray. Simply fold the count sheet to fit inside of the medical grade paper pouch and adhere to the side of your wrapped tray.



PBR-46 Count Sheet Holder

PBR-46

The PBR-46 Count Sheet Holder has a transparent poly front with a medical grade paper back. This product safely holds the count sheet on the outside of a wrapped tray. The poly front allows for quick visual confirmation that the count sheet is present.