

EndoPro-Cam

EPRO-001SK



Evaluate the Condition of Rigid Endoscopes with the Endo-Pro Cam®

The EndoPro-Cam® is a portable endoscopic video system that assists healthcare workers with inspection of their facilities' rigid endoscopes. The lightweight, ergonomic EndoPro-Cam® allows for image capture, recording and adapts to C-Mount couplers.

The EndoPro-Cam® comes with a camera stand, a Mini HDMI cable for secondary monitor usage, and a test pattern card, which allows the user to evaluate the resolution and color accuracy.

Specifications	
Display	4.3 inch LCD Monitor
Resolution	1080p High Resolution
Power	2 Rechargeable Li-ion batteries & Power Bank
Light Source	10W LED
Memory	32GB SD Card
Recording Format	MOV
Still Image Format	JPG



Easy-to-Use



Capture Images



Record Video



Lightweight, Ergonomic Unit



Rechargeable Batteries




4.3 Inch LCD Monitor


Meet AAMI Guidelines for Visual Inspection

AAMI ST79 (2017) Amendment 2: 8.2 Inspection of Instruments

Each time a medical device is processed, it should be visually inspected for cleanliness and integrity.

Enhanced inspection with magnification, borescopes, or other inspection methods to verify cleanliness and integrity may be used.

 hmark.com

 800.521.6224

 healthmark@hmark.com