


# Tear & Moisture Management Solutions

Healthmark offers a variety of solutions to help your department ensure trays are protected from damage and wet pack. Solutions include liners for inside of trays, under trays, autoclave shelves, heavy sets and low temperature sterilization. Interested in learning about these and other solutions? Contact us!

**hmark.com**  
**800.521.6224**



## For Inside or Under Trays

### UnderGuard™ Tray Liners


Line the inside of trays to protect instrumentation or use externally for tray corner protection and moisture-wicking.

### Dual UnderGuard™ Dry Mats

These moisture-wicking liners feature a double layer of tear resistance for optimum protection of instruments from dust, moisture and handling.

### UnderGuard™ Dry Advanced

Specially designed for moisture wicking and to prevent tears in wraps as well as to protect instruments internally. Surprisingly soft, but extremely durable and non-linting.



## For Inside of Trays



### Form-Fitting UnderGuard™ Tray Liners

Pre-cut tray liners that easily fold to form-fit into tray corners for protection of instrumentation.




## For Autoclave Shelves



### UnderGuard™ for Autoclave Shelves

Sized for use on autoclave shelves, this heavy duty shelf liner protects wraps from tears or damage and is moisture-wicking.




## For Under Trays



### OrthoShield™ Soaker Sheets

These moisture-wicking Soaker Sheets are specially developed for use under heavier orthopedic trays to protect wraps from sharp corners and uneven surfaces.



## For Low Temperature



### UnderGuard™ for VH202 Sterilization

Compatible with hydrogen peroxide sterilization. Use under the tray to protect wrap from sharp corners and uneven surfaces.