

Brand Name of Product	UnderGuard™ Dry Advanced
Generic Name of Product	Dry Mats
Product Code Number(s)	DA-2030PP, DA-2030GN, DA-1520PP, DA-1520GN, DA-1224PP, DA-1224GN, DA-1212PP, DA-1020PP, DA-1020GN, DA-1010PP, DA-1010GN
Intended Use	<ul style="list-style-type: none"> For use inside or outside of trays within a sterile barrier. Protects items while simultaneously providing moisture wicking against the dispersion of condensate.
Range of Applications for Product	Provides protection useful for trays and baskets.
Key Specifications of Product	<ul style="list-style-type: none"> Polypropylene (PP) Single use Non-sterile Non-linting.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	<ul style="list-style-type: none"> Temperature: 10- to 35 °C (50- to 95 °F) Humidity: 20–70% rh Store in dry areas away from sunlight. Not temperature controlled.
Packaging Contents	N/A
Shelf Life	Five (5) years from manufacturer date.

Instructions for Using Product	
Description of Use(s)	Protects items while simultaneously providing moisture wicking against the dispersion of condensate.
Preparation for Use	N/A
Diagrams (drawings, pictures)	N/A
Steps for Use of Product	Liners may be placed either under the tray or inside the bottom of the tray.
Interpretation of Test Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A
Special Warnings and Cautions	N/A
Disposal	<ul style="list-style-type: none"> Single use. Dispose of according to facility's policy.

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	<ul style="list-style-type: none"> Steam sterilization at: <ul style="list-style-type: none"> 132 °C (270 °F) for 4 minutes 135 °C (275 °F) for 3 minutes. Gravity steam sterilization at 121 °C (250 °F) for 30 minutes. Maximum tested time/temperature: 135 °C (275 °F) for 30 minutes.

	<ul style="list-style-type: none"> • Materially compatible with VH202 and EO sterilization.
Storage	<ul style="list-style-type: none"> • Temperature: 10- to 35 °C (50- to 95 °F) • Humidity: 20–70% rh • Store in dry areas away from sunlight. • Not temperature controlled.
Additional Information	N/A
Related Healthmark Products	Dual UnderGuard™ Dry Mats
Other Product Support Documents	Sterilization Brochure
Reference Documents	N/A
Customer Service Contact	Healthmark, A Getinge company 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com