

## Steriking® Reinforced Rolls

Developed for instruments requiring stronger packaging.



### The smart choice for heavy instruments and trays.

- Steriking® reinforced rolls are perfectly suited to demanding applications such as orthopedic instruments, tray sets and stainless steel holloware.

# Steriking® Reinforced Rolls



Steriking® reinforced rolls allow optimum choice and flexibility for processing of each demanding application.

- Strong materials keep the instruments sterile during transportation and storage
- Long shelf life after sterilization
- Savings in material consumption and less waste compared to tray wrap



## Steriking® Reinforced Rolls

Steriking Code	Width of the roll (inches)	Examples of packing applications. This is not an all encompassing list.
Heat-Sealable pouches:		
RR41	4 inches	Bone elevators, bone curettes, bone hooks, Cobb elevator, orthopedic screwdrivers, laparoscopic instruments, rake retractors, small mallets, chondro picks, trocars, osteotomes, spinal curettes, dental instrument cassettes, ash dental forceps
RR43	8 inches	Shoulder retractors, hip retractors, kidney bowls, bone breakers, bone nibblers, hammers, coeloscopy forceps, bone curettes, spider beach retractor, ortho graspers, spinal Kerrison rongeurs, medium ortho mallets, angled gelpi, vaginal speculums, Wrigley forceps, Neville Barnes forceps, laparoscopic needle holders, laminar spreaders, weight and chains, implant trays
RR45	12 inches	Orthopedic heavy mallets, stainless steel jugs, light handles, large Charnley hip curette, pin holders and pins, micro-instrument sets, orthopedic, small basins, endoscope baskets, anaesthetic trays
RR47	16 inches	PP instrument baskets, PP thermoform trays, stainless steel jug, large light handles, large kidney dishes, foot support, instrument sets, round prep basins, retractors, Deaver retractors, Charnley bow retractors
RR49	20 inches	Instrument sets, endo eyes, rinse bowls, basins





# Steriking® Reinforced Rolls



Steriking® reinforced rolls allow optimum choice and flexibility for processing of each demanding application.




Steriking® reinforced rolls are the optimum solution for difficult applications and provide an innovative solution to common issues faced daily in sterile processing departments.

High performance demands high standards, increase your department performance with Steriking sterile barriers.



# Steriking® Reinforced Rolls

Steriking® renowned quality



## High protection



- Strong medical paper and Steriking's Multi-X9 film for superior protection against punctures and tearing
- Outstanding microbial barrier and protection against contamination

## Efficient packaging

- Using the Steriking® reinforced rolls allow for 2-5 times faster packaging compared to tray wrap or container
- Sealing temperature will be the same as with the other STERIKING paper/laminate products
- Transparent Multi-X9 film allows to see the packed content

## Aseptic presentation

- Clean fiber-free opening



Instrument symbol marking on the roll allows users to define packing the instruments right way to aseptically present.



Wipak Oy  
Wipaktie 2  
FI-15560 Nastola

At Wipak, we take Sustainability seriously

- using green energy in production, recycling our production waste, using shipping boxes that are made from recycled materials, making products and innovations that have long shelf life before and after sterilization.



*We care that you pack safely!*  
**www.steriking.com e-mail: steriking@wipak.com**