

Getinge Assured

Mini Bowie-Dick Test Pack

The Mini Bowie-Dick Test Pack is a single-use device designed to detect the presence of residual air in pre-vacuum steam sterilizers operating at 132°C and 134°C for 3.5 minutes. The indicator card within the Mini Bowie-Dick Test Pack will demonstrate a uniform color change from purple to green when proper sterilization conditions are met and no air is present.



Performance Characteristics

The indicator card within the Mini Bowie-Dick Test Pack will demonstrate a uniform color change from purple to green when proper sterilization conditions are met and no air is present. If enough air is present to create a 2°C (+1°/-0°C) temperature difference between the center of the towel pack, as identified in ANSI/AAMI/ISO 11140-5 and the drain temperature at the beginning of the final one minute of a three and half minute cycle the indicator card within the Mini Bowie-Dick Test Pack will demonstrate a non-uniform color change.

Key Features

- Complete color change - 3.5 minutes at 132°C and 134°C
- Less paper in each pack which reduces the amount of paper waste
- Test pack is 50% smaller than other Bowie-Dick test packs

Physical Characteristics

- **Lamination:** Top and bottom
- **Dimensions Card:** 3.97" x 2.47"
- **Packaging:** 30 test packs

Storage

The storage conditions are printed on the product packaging

Shelf Life

24 months. The expiry date is printed on the product packaging

Standards Classification

FDA Cleared Air Removal Indicator
November 26, 2013, K132372

Quality Statement

Confidence in Getinge is the most important quality criteria. This is the hallmark of all our external and internal commitments, activities and products. Products and services supplied by Getinge conform to the agreed terms and expectations. The achievement of these quality goals is the basis for continued competitive and successful enterprise.

Ordering Information

Getinge Assured

Mini Bowie-Dick Test Pack

Order Number

61301606667

Getinge is a registered trademark of Getinge AB or its subsidiaries or affiliates in the United States or other countries. • Copyright 2019 Getinge AB or its subsidiaries or affiliates. • All rights reserved. • MCV00054059 REV B

Manufactured by: SteriTec Products, Inc. 74 Inverness Drive East, Englewood, CO 80112 USA
Sales Office: Getinge USA Sales, LLC 1777 East Henrietta Road, Rochester, NY 14623 USA

www.getinge.com