

Equipment Name		Model	
----------------	--	-------	--

SonoCheck Location*	Result	Water Temperature	Program Cycle	Comments	Tester Initials	Date

Testing of the ultrasonic should be done whenever the equipment is used.

1	2	3
4	5	6
7	8	9

\*The image on the left represents a standard basin, divided into 9 zones. Please refer to the grid to indicate the SonoCheck Location (1–9).



**FAIL**

Blue/Green



**FAIL**

Light Blue/Green



**PASS**

Yellow



**PASS**

Clear (rare)