

AAMI LEVEL 4 Decontamination Gown

Comfort without compromise

The AAMI Level 4 Decontamination Gown is a full coverage gown designed specifically for the demands of sterile processing workflows, where exposure to fluids and strike-through risks are a daily concern.

Compliant with **ANSI/AAMI PB70:2022 Level 4** fluid barrier requirements, it helps protect against fluid and microorganism exposure during device sorting and decontamination activities. The gown features a reinforced diamond front panel designed to provide enhanced coverage and strike-through resistance in high-exposure areas.

- ✓ **AAMI Level 4** fluid barrier protection
- ✓ **Purpose-built** for the demands of decontam
- ✓ **Reinforced front panel** for enhanced coverage
- ✓ **40% lighter** construction for wearability*
- ✓ Designed for **easier donning & doffing**



Reinforced Front Panel

Diamond-shaped coverage designed to help protect key exposure areas



Secure Fit Features

Hook-and-loop neck closure, thumb loops, and pull tab support a controlled fit & removal



Adjustable Waist Closure

Waist ties provide a secure fit that adjusts to you and stays in place