

Deliver a Clear View with **EndoPro-Cam[®]**

Inspect and document rigid scope
image quality for maintenance and
quality assurance activities

The EndoPro-Cam helps bring greater clarity to rigid scope inspection. This portable 1080p HD imaging system enables users to visually assess rigid endoscope image quality and capture images or video to support documentation.

Complete with a standardized test pattern card for visual assessment using a consistent reference image, C-Mount compatibility, and a lightweight, ergonomic design, the EndoPro-Cam supports consistent, convenient inspection during maintenance and quality assurance activities.

hmark.com
800.521.6224



Identify Image Quality Issues

Visually assess rigid endoscope image output for visual defects and to document observations.

Document with Confidence

Capture still images and video to support consistent documentation and communication.

Inspect with Ease

Lightweight, ergonomic, and rechargeable for convenient inspection throughout the workday.

Assess Image Quality

Use the included standardized test pattern card as a consistent visual reference when observing image resolution, color reproduction, and other image characteristics.



Learn More

Scan to explore how the EndoPro-Cam can support your rigid scope inspection process.